



Guide to PROTEIN EXPRESSION

Index

| Media Preparation | 4 | |
|------------------------------|--------------------|--|
| Overnight culture prepar | ation 8 | |
| Overexpression | 13 | |
| Cell Harvesting and Lysi | ng 17 | |
| Centrifugation and Purifi | cation 23 | |
| Conclusion | 28 | |
| Ordering Information | 30 | |
| CAPPRondo Hotplate Stirrer | 30 | |
| CAPPSolo Single Channel P | Pipette 30 | |
| BluCAPP Centrifuge Tubes | 31 | |
| Expellplus Low Retention Fil | ter Tips 31 | |
| CAPPRondo Platform Shake | er 31 | |
| CAPPRondo Microcentrifuge | e 32 | |
| CAPPRondo Magnetic Hotpl | late Stirrer 32 | |
| Expell Cryoboxes - plastic | 32 | |
| Expell Cryoboxes - cardboar | d 32 | |
| CAPPController | 32 | |
| CAPPHarmony Serological F | Pipetttes 33 | |
| Expel / Expel Secure Microc | entrifuge tubes 33 | |



Protein expression and harvesting has brought about some exciting discoveries in life sciences.

With the vast amount of information now available for us to extract from the proteome, these techniques unlock a wealth of knowledge that can benefit the future of humanity and the natural world.

CAPP supplies the innovative tools that give scientists an edge in this protein expression gold rush. By empowering scientists to do their best work, CAPP products improve humanity's understanding and encourage careful use of the vast ocean of proteins that are available to study.

In this brief guide, we will outline how CAPP products can be leveraged in protein expression and harvesting workflows by demonstrating their applicability to a simulated protein expression workflow^{*}.

All state-of-the-art laboratory products featured in this guide are the result of CAPP's 20+ years of experience in manufacturing world-class laboratory instruments.

To learn more about the different protein expression and harvesting solutions from CAPP, you can visit our website (<u>www.capp.dk</u>) and select the product of interest from the navigation menu.

*For detailed protein expression protocols please refer to published scientific literature.

Media Preparation

Good media preparation is the underlying foundation for successful protein expression.

In this section, we outline how the 23-litre CAPPRondo Hotplate Stirrer gives you the edge in your media preparations.



Rehydrate the LB broth by adding an appropriate amount of solvent, according to the manufacturer's instructions.

By using manufacturer-recommended solvents, you lower the chances of poorly prepared media, which may impact negatively on cell growth and result in lower protein yields.

Helpful tips:

- Only use high-purity solvents without physical and/or chemical impurities
- Always use clean containers for your media preparation
- Ensure all glassware and plasticware is fully rinsed and detergent-free

Results-oriented Innovation with CAPP

Using the CAPPRondo Hotplate Stirrer, thoroughly dissolve the media suspension.

The CAPPRondo Hotplate Stirrer provides even heating for quick dissolution of the media suspension. With a fully adjustable timer and variable speed settings, the CAPPRondo Hotplate Stirrer can be pre-programmed to run continuously with minimal operator supervision.

CSPP Rondo

HOTPLATE STIRRER

Features:

| Speed Range | 200-2.200 RPM in steps of 10 RPM |
|----------------------|----------------------------------|
| Temperature Range | From ambient to 550°C |
| Load Capacity | 20L |
| Motor Type | Brushless DC motor |
| Top Plate Material | Nano-crystalline glass ceramic |
| Top Plate Dimensions | 180x180mm |
| Time Setting | 1-999 minutes or infinite |
| Programmable Modes | Up to 99 individual protocols |
| Pulse Mode | Programmable, from 30 to 99 secs |
| Dimensions (WxDxH) | 220x330x115mm |
| Weight | 4.5kg |
| Protection Class | IP 21 |
| | |

Once the LB broth is fully rehydrated, sterilize it by autoclaving at 121°C for 15 minutes.

After autoclaving, pour the sterile media into sterile containers in a laminar-flow cabinet. Sterilization helps to lower the chances of contaminating the culture with potentially unwanted microorganisms.

Helpful tips:

- Use a clean pair of gloves to handle the newly sterilized media and change gloves regularly
- Always ensure that the laminar-flow cabinet is serviced regularly by qualified personnel
- Remember to use appropriate tools for handling the hot media containers/flasks to protect your fingers from burning



Overnight culture preparation

Optimize your overnight culture prep with CAPP's superior quality instruments and consumables.

Transfer the desired volume of the LB broth into sterile 50mL BluCAPP Centrifuge Tubes.

BluCAPP Centrifuge Tubes are preferable to glass flasks because of their ease of handling and excellent compatibility with shakers and vortex mixers. With easy-to-read graduations and ³/₄-turn closures, CAPP's BluCAPP Centrifuge Tubes help scientists to speed up their protein expression workflow.



BluCAPP CENTRIFUGE TUBES

Features:

The polypropylene-constructed tubes are autoclavable and tested to endure up to 17,000 RCF

Ridged leak-proof caps provide an optimal grip which, when combined with a ³/₄-turn closure, allows single-handed operation and significantly increases speed

Using an autoclavable CAPPSolo Single Channel Pipette (20-200 µL) with presterilized CAPP ExpellPlus Low Retention Tips, transfer an adequate amount of antibiotic into a BluCAPP Centrifuge Tube to suppress the growth of unwanted microorganisms.

CAPPSolo Single Channel Pipettes provide dependable pipetting accuracy and precision, allowing you to transfer the exact volume of antibiotic every time.



SINGLE CHANNEL PIPETTES

Features:

Exchangeable colour inserts allow easy pipette identification, either by user or volume capacity

Fine-tuned design ensures high accuracy and long-term calibration retention

Light pipette chassis reduces any risk of repetitive strain injuries

Counter-mechanism prevents accidental rotation while pipetting

Withstands the toughest laboratory conditions including chemical exposure, drops and autoclaving

Tip cone allows universal tip compatibility

With the CAPPSolo Single Channel Pipette and Expellplus Low Retention Tips, inoculate the antibiotic-treated media in 50mL BluCAPP Centrifuge Tubes.

This is done by transferring *E. coli* cells expressing the chosen protein from a stock culture. The comfortable, lightweight body of CAPPSolo Single Channel Pipettes makes it easy to keep your hands perfectly steady as you transfer the inoculum to BluCAPP Centrifuge Tubes.



With CAPP you culture with confidence

Mix the culture with the CAPPRondo Platform Shaker and incubate overnight at 37°C for the cells to grow and express the protein of interest.

This excellent shaker has multiple attachments, making it universally compatible with many types of culture containers. The CAPP Platform Shaker includes an auto-restore feature that enables the previous programme to continue, even in the event of a power outage.



PLATFORM SHAKER

Features:

| Shaking Movement | Orbital (Ø20mm) / Linear |
|-----------------------|---|
| Max. Loading Capacity | 7.5kg |
| Speed Range | 50-300 RPM |
| Motor Type | Brushless DC motor |
| Time Settings | 0 to 99 hours and 59 minutes |
| Dimensions (WxDxH) | 408x350x99mm |
| Platform Size | 330x330mm |
| Protection Class | IP 21 |
| Programmable Features | Up to 6 individual protocols, each with |
| | additional 9 sub-programmes |
| Programmable Looping | YES |
| RS 232 Interface | Optional |
| | |



Overexpression

Prepare for successful overexpression the CAPP way.

The following day, inoculate the main culture with the overnight culture and antibiotic using an autoclavable CAPP pipette with 1mL sterilized tips.

CAPP autoclavable pipettes are easy to sterilize in an autoclave. These, when combined with 1mL CAPP sterilized tips manufactured in a cleanroom facility, will minimize the chances of contaminating the culture.

CAPP Pipette Tips are constructed from premium-grade polypropylene and have industry-leading visual clarity.



Universal tip compatibility

Compatible with a large number of pipette brands that are available on the market.

High quality filter

Made of HDPE without cellulose additives, these filters are inert and feature a pore size that is 3 to 4 times smaller than several other filter brands. Researchers working with protein expression samples can have peace of mind knowing that they have a reliable filter to prevent cross contamination.

Certified DNA, RNA & Pyrogen-free

The ExpellPlus filter tips are tested to conform to strict bio-load specifications. Batch certificates are available on request.

Low retention polymer

The ExpellPlus low retention polymer virtually eliminates any liquid retention on the tip, ensuring greater repeatability of liquid transfers and less wastage of expensive reagents. Using the CAPPRondo Platform Shaker, agitate the newly cultured medium at the desired temperature in a climate-controlled environment such as an incubator or hood.

The CAPPRondo Platform shaker is designed to make small-scale protein expression workflows much easier.





PLATFORM SHAKER



After completing overexpression, transfer the culture into 1.5mL CAPP Cryo Tubes in preparation for SDS-PAGE.

CAPP Expell 1.5mL Cryo Tubes are available in an exclusive, secure design that minimizes any risk of accidental opening or sample loss.

Lip seal cap design ensures that the tubes are leak-proof

QuickSeal[™] design facilitates single-handed operation. Tighten or open it with only half a turn!

SuperStraight[™] design fits into both 9x9 and 10x10 cryoboxes. No surface transition between the tube and the cap for minimal risk of contamination during storage

1.0 0.5 H Coloured cap inserts Available in 5 colours. Made of polypropylene for easy tube recognition

External thread SuperStraight™ design gives up to 23% extra storage capacity when used with the 10x10 cryoboxes

High quality graduation marks and writing area for easy sample identification

> Star-foot bottom for more stability and better operation when stored in selected cryoboxes

EXPell Cryo Tubes

Cell harvesting and lysing

Harvest more cells with CAPP instruments and consumables.

Using the CAPPRondo Microcentrifuge, spin down the CAPP Cryo Tubes into separate pellet from the supernatant.

The elegantly-designed CAPPRondo Microcentrifuge features a powerful, ultra-quiet motor that makes quick work of spinning down tubes with induced cells. With its small footprint, the CAPPRondo Microcentrifuge is perfect for space-limited settings.



CSPP Rondo

MICROCENTRIFUGE

Features:

CR-68

| Motor Type | DC motor |
|--------------------|---|
| Rotor Capacity | 8 x 1.5/2.0mL & 2-place PCR strip rotor |
| Speed | 2000-6000 RPM |
| | 2000 x g (MAX RCF) |
| Time | No |
| Dimensions (WxDxH) | 162 x 157 x 116 mm |
| Weight | 1.1kg |
| | |

CR-68X

| Motor Type | Brushless DC motor |
|--------------------|---|
| Rotor Capacity | 8 x 1.5/2.0mL & 2-place PCR strip rotor |
| Omenad | 2000-6000 RPM |
| Speed | 2000 x g (MAX RCF) |
| Speed Accuracy | ± 25 RPM |
| Time | 0 to 25 mins & infinite |
| Dimensions (WxDxH) | 162 x 157 x 116 mm |
| Weight | 1.1kg |
| - | |

After centrifuging the samples, carefully remove the supernatant from the Cryo Tube using a CAPPSolo Single Channel Pipette and CAPP Tips.

CAPPSolo Single Channel Pipettes feature an easy-to-hold design that gives users perfect control of their pipetting experience, preventing disturbance of the pellet.



SINGLE CHANNEL PIPETTES

Features:

Made in Germany

Fully autoclavable

Robust, yet extremely lightweight

Ergonomically comfortable with curved finger rest and soft handle

Exchangeable coloured drop inserts for easy identification

Low cost, quality design with the proven legacy of CAPP metallic shafts!

Extremely reliable tip fit/eject

Add an appropriate volume of fresh media to the CAPP Cryo Tube to resuspend the cell pellet.

Using your CAPPSolo Single Channel Pipette and a fresh tip, gently mix the sample by slowly pipetting up and down three times.

For even more efficient processing, there are colour-coded cap inserts for the CAPP Expell Cryo Tubes that make it easier to categorize samples at various stages of the processing workflow. Coming in red, blue, green, white, and yellow these cap inserts facilitate easy sample sorting and also add a touch of colour to your workbench.



Features:

One of the only external-thread tubes on the market that is suitable for storage in 10x10 boxes, providing up to 23% extra storage capacity

Sturdy star-foot bottom enables opening and closing of the tube without taking it out of the box

Suitable for storage in 9x9 and 10x10 boxes

Available in 5 sizes: 0.5mL, 1mL, 2mL, 4mL and 5mL

Lip seal cap ensures it is 100% leak-proof without using O-rings

QuickSeal design for single-handed operation, enabling you to tighten or open the cube with only half a turn

Exchangeable colour-coded caps for easy identification in 5 colours



Transfer the resuspended induced cells to 50mL CAPP sterilized Centrifuge Tubes and store them at -80°C.

For ease of storage, we recommend packing your CAPP Centrifuge Tubes in CAPP Cryoboxes as they make organizing your samples quick and tidy. Although universally compatible with many other brands of cryo tubes, when used with CAPP's Cryo Tubes they provide further unique, space-saving advantages due to our unique design features.

Features:

Unique cryoboxes line with 9x9 and 10x10 dividers

Made of a high-quality durable polypropylene which is suitable for storing at -90°C and autoclaving at +121°C for 20 minutes

Star-foot lock at the bottom of each well enables single-handed opening of the tube without taking it out of the box

5 sizes of EXPell Cryoboxes: height: 32, 40, 50, 80 and 90mm perfectly fits all sizes of EXPell Cryotubes

Numerical coding on both the lid and the bottom of the box for easy identification and recording

Universally compatible

EXPell Cyoloxes PLASTIC CRYOBOXES

If you are looking for an environmentally friendly option for your cryo storage, CAPP's Cardboard Cryoboxes are a great alternative, due to being fully recyclable after use.

CAPP Expell Cardboard Cryoboxes are the ultimate solution for environmentally friendly cryogenic storage. These highly flexible storage units provide exceptional value to both small life science labs and extensive research centers that are heavily invested in cryogenic storage of samples and reagents.

Features:

Constructed from high-quality cardboard with a water-repellent finish for extra durability

Thermostable, recyclable cardboard

Removable grid allows more flexibility for your storage needs

Supplied in 9x9 and 10x10 grid configuration options

Write-on surface makes identifying your samples a lot easier

Chemically-resistant to alcohols and mild organic solvents

EXPell Constant Cardboard Cryoboxes

Centrifugation and Purification

Optimize your protein harvesting for success by choosing the best tools from CAPP's toolbox.

With the aid of the ceramic-coated CAPP Magnetic Hotplate Stirrer, prepare lysis buffer(s).

This programmable hotplate stirrer has intuitive controls that make it easy to adjust key parameters such as speed and time. With the last-run memory feature, it is easy to continue with the last-used programme; therefore minimizing any need to fiddle with settings every time.

With a 10L capacity, the CAPP Magnetic Hotplate Stirrer handles a massive load for device of its size.





MAGNETIC HOTPLATE STIRRER

Features:

| Max. Load Capacity | 10L |
|-------------------------|---|
| Speed Setting | 200 - 2.200 RPM, in steps of 10 RPM |
| Temp. Range | from Ambient to 320°C |
| Time Setting | 1-999 mins and continuous |
| Digital Display | Yes. Shows Speed/Time/Actual/Set Temp. alternately |
| Safe Temp. Circuit | Yes |
| Plate Diameter | 140 mm |
| Weight | 1.4kg |
| Dimensions (WxDxH) | 248x156x104mm |
| Protection Class | IP 21 |
| Accessories included | CRS-21H-P and CRS-21H-S |
| Pulse Mode | Run-Pause-Run. Interval of 30 seconds |

Using the advanced CAPP Pipette Controller and a CAPP sterilized Serological Pipette, transfer the appropriate amount of lysis buffer into a conical tube with induced cells.

At just 167g, the CAPP Pipette Controller is a lightweight, ergonomically pleasing device that enables rapid, fatigue-free, high-accuracy liquid transfer (±2%).

Features:

Powerful motor can fill a 50mL pipette in less than 10 seconds

Gravity delivery and blow-out dispensing mode selection

Fully adjustable speed-setting

Exchangeable 0.2 µm PTFE disc filter prevents overflow

Vapour outlet for venting out liquids or fumes

Rechargable Lithium battery up to 12 hours of continuous use

Fully autoclavable nose cone

UV-resistant housing

Complies with CE and RoHS



PIPETTE CONTROLLER

CappAid"

ddk



When the CAPP Pipette Controller is used together with cotton-plugged, pyrogen-free CAPP sterilized Serological Pipettes, scientists can minimize any risk of sample contamination.

Designed for accurate, rapid fluid dispensing, these transparent, medical-grade polystyrene pipettes are sterile, pyrogen-free and cotton-plugged for many liquid handling applications.

Features:

Sharp, highly legible permanent graduations prevent ambiguity and guarantee precise reading and dispensing to +1-2%

Each pipette comes with an effective filter barrier to prevent any overflow of liquid into the pipette

Reverse graduations show the volume removed or remaining on all pipettes larger than 1mL

The colour-coded mouthpiece fits all major pipettes

After lysis, separate the soluble and insoluble fractions by centrifugation using CAPPRondo Microcentrifuge Tubes.

With the aid of the 200 microliter CAPPSolo single-channel micropipette, separate the supernatant and pellet into separate 1.5mL CAPP Microcentrifuge Tubes for storage, pending further processing and purification.

CAPP's 1.5mL Microcentrifuge Tubes are available in standard (Expell) and low-retention (ExpellPlus) versions.



Black tubes For light-sensitive applications **Single handed** open/close operation **Unique Click-Lock lid** Leak-proof, minimizes **Perfect sealing and** the risk of opening evaporation protection during operation Temperature range from -80 to 121°C. **Polished interior Coloured tubes** ensures good sample For easy sample visibility differentiation

Conclusion

While the simulated protein expression workflow in this guide may differ from real world scenarios, it demonstrates that successful protein expression and harvesting benefit from a clear understanding of the scientific principles underlying the procedures, as well as good technical skill for execution.

With CAPP products, every scientist can get excellent results from their protein expression workflows, regardless of their skill level. Every product highlighted in this short guide is designed to make science easy, fun, and reliable by offering intuitive operation, outstanding quality, and long-lasting performance.



With CAPP products, every scientist can get excellent results from their protein expression workflows, regardless of their skill level

39

Ordering Information

CIPP Rondo HOTPLATE STIRRER



| Cat. No. | Description |
|-----------------|---|
| CRS-22H | CAPPRondo Hotplate Stirrer 2.200 RPM, 550°C ceramic, 20L capacity, 220V |
| CRS- 22H-110 | CAPPRondo Hotplate Stirrer 2.200 RPM, 550°C ceramic, 20L capacity, 110V |
| CRS-22H-P | P1000 temperature probe for CAPPRondo Hotplate Stirrer |
| CRS-22H-S | Stand for P1000 temperature probe CAPPRondo Hotplate Stirrer |
| | |

CSPP Solo

SINGLE CHANNEL PIPETTES

CAPP Solo Fixed Volume

| Cat. No. | Description |
|-----------|-------------------------------------|
| S5-1F | CAPPSolo pipette, fixed vol. 5µL |
| S10-1F | CAPPSolo pipette, fixed vol. 10µL |
| S20-1F | CAPPSolo pipette, fixed vol. 20µL |
| S25-1F | CAPPSolo pipette, fixed vol. 25µL |
| S50-1F | CAPPSolo pipette, fixed vol. 50µL |
| S100-1F | CAPPSolo pipette, fixed vol. 100µL |
| S200-1F | CAPPSolo pipette, fixed vol. 200µL |
| S250-1F | CAPPSolo pipette, fixed vol. 250µL |
| S500-1F | CAPPSolo pipette, fixed vol. 500µL |
| S1000-1F | CAPPSolo pipette, fixed vol. 1000µL |
| S2000-1F | CAPPSolo pipette, fixed vol. 2mL |
| S2500-1F | CAPPSolo pipette, fixed vol. 2,5mL |
| S5000-1F | CAPPSolo pipette, fixed vol. 5mL |
| S10000-1F | CAPPSolo pipette, fixed vol. 10mL |

......

CAPP Solo Variable Volume

CAPP Solo StarterKit

| Cat. No. | Description |
|----------|--|
| S2-1 | CAPPSolo pipette, variable vol. $0.1 - 2.5 \mu L$ |
| S10-1 | CAPPSolo pipette, variable vol. $0.5 - 10 \mu L$ |
| S20-1 | CAPPSolo pipette, variable vol. 2 – 20 μ L, but XL tip fit |
| S50-1 | CAPPSolo pipette, variable vol. 5 – 50µL |
| S100-1 | CAPPSolo pipette, variable vol. 10 – 100µL |
| S200-1 | CAPPSolo pipette, variable vol. $20 - 200 \mu L$ |
| S1000-1 | CAPPSolo pipette, variable vol. 100 – 1000µL |
| S5000-1 | CAPPSolo pipette, variable vol. 0.5 – 5mL |
| S10000-1 | CAPPSolo pipette variable vol. 1 – 10mL |



| Cat. No. | Description |
|-----------|--|
| SK-03-S01 | 3 pipettes S10-1 /S100-1 /S1000-1, C-05 stand, 3 racks of tips and 1 timer |
| SK-03-S02 | 3 pipettes S20-1 /S200-1 /S1000-1, C-05 stand, 3 racks of tips and, 1 timer |



| Cat. No. | Description |
|----------|---|
| 5100014C | BluCAPP Centrifuge Tubes, 15mL, bulk, 20x25 pcs. |
| 5100015C | BluCAPP Centrifuge Tubes, 15mL, pre-sterile, 20x25 pcs. |
| 5100019C | BluCAPP Centrifuge Tubes, 15mL, cardboard rack, pre-sterile, 20x25 pcs. |
| 5100049C | BluCAPP Centrifuge Tubes, 50mL, bulk, 20x25 pcs. |
| 5100050C | BluCAPP Centrifuge Tubes, 50mL, pre-sterile, 20x25 pcs. |
| 5100059C | BluCAPP Centrifuge Tubes, 50mL, cardboard rack, pre-sterile, 20x25 pcs. |
| 5100060C | BluCAPP Centrifuge Tubes, 50mL, pre-sterile, self- standing, 20x25pcs. |
| | |



LOW RETENTION FILTER TIPS

| Cat. No. | Description |
|----------|--|
| 5030030C | ExpellPlus 10µL, pre-sterile w/ filter, hinged racks, 5x10x96 pcs. |
| 5030060C | ExpellPlus 10µL XL, pre-sterile w/ filter, hinged racks, 5x10x96 pcs. |
| 5030061C | ExpellPlus 10/20µL XL, extra narrow, pre-sterile w/ filter, hinged racks, 5x10x96 pcs. |
| 5030062C | ExpellPlus 20µL, pre-sterile w/ filter, hinged racks, 5x10x96 pcs. |
| 5030066C | ExpellPlus 100µL, pre-sterile w/ filter, hinged racks, 5x10x96 pcs. |
| 5030090C | ExpellPlus 200µL, pre-sterile w/ filter, hinged racks, 5x10x96 pcs. |
| 5030115C | ExpellPlus 300µL, pre-sterile w/ filter, hinged racks, 5x10x96 pcs. |
| 5030150C | ExpellPlus 1000µL, pre-sterile w/ filter, hinged racks, 4x8x96 pcs. |
| | |





| Cat. No. | Description |
|-------------|--|
| CRP-3X | CAPPRondo Platform Shaker ø 20mm, with linear and orbital motion, 7.5kg capacity |
| CRP-U | Universal Platform with Adjustable Bars (Flask capacity 9 x 250mL, 6 x 500mL) |
| CRP-D | Disk Attachment Base Plate with Slip- resistant tray |
| CRP-C9-250 | Base Plate with 9 Plastic Clamps: 9 x 250mL |
| CRP-C9-500 | Base Plate with 9 Plastic Clamps: 9 x 500mL |
| CRP-F | Funnel Attachment |
| CRP-CS9-250 | Base Plate with 9 Stainless Steel Clamps: 9 x 250mL |
| CRP-CS9-500 | Base Plate with 9 Stainless Steel Clamps: 9 x 500mL |
| CRP-CS-100 | Stainless Steel Clamp for 100mL Flask |
| CRP-CS-150 | Stainless Steel Clamp for 150mL Flask |
| CRP-CS-250 | Stainless Steel Clamp for 250mL Flask |
| CRP-CS-500 | Stainless Steel Clamp for 500mL Flask |
| CRP-CS-1000 | Stainless Steel Clamp for 1000mL Flask |

Note: CAPP offers a wide portfolio of high-quality tips, i.e. Bulk Filter Tips.

See all options under: https://www.capp.dk/categories/pipette-tips



Cat. No. Description

| CR-68 | CAPPRondo Microcentrifuge 6.000 RPM/2.000g |
|--------|--|
| CR-68X | CAPPRondo Microcentrifuge w/ adjustable speed and timer function, max. 6.000 RPM/2.000g |



MAGNETIC HOTPLATE STIRRER

| Cat. No. | Description |
|-------------|--|
| CRS-22H | CAPPRondo Hotplate Stirrer 2.200 RPM, 550°C ceramic, 20L capacity, 220V |
| CRS-22H-110 | CAPPRondo Hotplate Stirrer 2.200 RPM, 550°C ceramic, 20L capacity, 110V |
| CRS-22H-P | P1000 temperature probe for CAPPRondo Hotplate Stirrer |
| CRS-22H-S | Stand for P1000 temperature probe CAPPRondo Hotplate Stirrer |

EXPell Cryoboxes PLASTIC CRYOBOXES

| Cat. No. | Description |
|----------|---|
| 5040510C | Expell Plastic Cryobox , 10x10 cell grid for 0.5mL cryovials, case of 36 pcs. |
| 5041010C | Expell Plastic Cryobox, 10x10 cell grid for 1,0mL cryovials, case of 36 pcs. |
| 5042010C | Expell Plastic Cryobox, 10x10 cell grid for 2,0mL cryovials, case of 36 pcs. |
| 5044010C | Expell Plastic Cryobox, 10x10 cell grid for 4,0mL cryovials, case of 36 pcs. |
| 5045010C | Expell Plastic Cryobox, 10x10 cell grid for 5,0mL cryovials, case of 36 pcs. |
| 5042022C | Expell Plastic Cryobox, 9x9 cell grid for 1,0 - 2,0mL cryovials, 130x130x50 mm, case with 36 pcs. |
| 5044020C | Expell Plastic Cryobox, 9x9 cell grid for 4,0mL cryovials, 130x130x80 mm, case with 36 pcs. |
| 5045020C | Expell Plastic Cryobox, 9x9 cell grid for 5,0mL cryovials, 130x130x90 mm, case with 36 pcs. |



CARDBOARD CRYOBOXES

| Cat. No. | Description |
|----------|---|
| 5042030C | Cardboard Cryobox B9, 133x133x50 mm, 9x9 cell grid, case with 36 pcs. |
| 5042020C | Cardboard Cryobox B10, 133x133x50 mm including 10x10 cell grid, case with 36 pcs. |





PIPETTE CONTROLLER

| Cat. No. | Description | |
|----------|--------------------------|--|
| Т100-В | CAPPTempo Blue | |
| T100-G | CAPPTempo Green | |
| T100-R | CAPPTempo Red | |
| T100-V | CAPPTempo Violet | |
| T100-CS | CAPPTempo Charging Stand | |



| Cat. No. | Description | Package | Colour |
|----------|---------------------------------------|--|--------|
| SP-1-C | CAPPHarmony Serological pipette 1mL | individually wrapped, 500 pcs. | • |
| SPB-1-C | CAPPHarmony Serological pipette 1mL | bulk, 25 pcs./bag, 40 bags/case, (1000 pcs.) | • |
| SP-2-C | CAPPHarmony Serological pipette 2mL | individually wrapped, 500 pcs. | |
| SPB-2-C | CAPPHarmony Serological pipette 2mL | bulk, 25 pcs./bag, 40 bags/case, (1000 pcs.) | |
| SP-5-C | CAPPHarmony Serological pipette 5mL | individually wrapped, 200 pcs. | |
| SPB-5-C | CAPPHarmony Serological pipette 5mL | bulk, 25 pcs./bag, 20 bags/case (500 pcs.) | |
| SP-10-C | CAPPHarmony Serological pipette 10mL | individually wrapped, 200 pcs. | • |
| SPB-10-C | CAPPHarmony Serological pipette 10mL | bulk, 25 pcs./bag, 20 bags/case (500 pcs.) | • |
| SP-25-C | CAPPHarmony Serological pipette 25mL | individually wrapped, 150 pcs. | • |
| SPB-25-C | CAPPHarmony Serological pipette 25mL | bulk, 25 pcs./bag, 8 bags/case (200 pcs.) | |
| SP-50-C | CAPPHarmony Serological pipette 50mL | individually wrapped, 100 pcs. | |
| SPB-50-C | CAPPHarmony Serological pipette 50mL | bulk, 25 pcs./bag, 4 bags/case (100 pcs.) | |
| SP-100-C | CAPPHarmony Serological pipette 100mL | individually wrapped , 50 pcs. | |



MICROCENTRIFUGE TUBES

| Cat. No. | Description |
|----------|--|
| 5100500C | Expell Microcentrifuge tubes 0.5mL, bulk, 10x1000 pcs. |
| 5001500C | ExpellPlus Microcentrifuge tubes 1.5mL, bulk, low retention, 20x500 pcs. |
| 5101500C | ExpellSecure Microcentrifuge tubes 1.5mL, bag, 20x500 pcs. |
| 5101501C | Expell Microcentrifuge tubes 1.5mL, bag, 20x500 pcs. |
| 5101590C | ExpellSecure Microcentrifuge tubes 1.5mL, bulk, assorted colours, 20x500 pcs. (magenta, blue, yellow, green, orange) |
| 5101550C | ExpellSecure Microcentrifuge tubes 1.5mL, bulk, black, 20x500 pcs. |
| 5102000C | ExpellSecure Microcentrifuge tubes 2.0mL, bag, 16x500 pcs. |
| 5102001C | Expell Microcentrifuge tubes 2.0mL, bag, 10x500 pcs. |

With CAPP, you have a partner committed to delivering outstanding scientific results every single time...

> CAPP, truly accurate, truly robust!



Uthleber Weg 14 99734 Nordhausen Germany

P: +49(0)3631/65242-0 F: +49(0)3631/65242-90 E: info@cappahn.com www.capp.dk



© AHN Biotechnologie GmbH 2021