

If you are looking for a simple, yet effective solution to your reagent storage applications, look no further, **CAPP Reagent Bottles** have got you covered. Drawing from our extensive manufacturing expertise, our Reagent Bottles tout a high-quality design that offers versatility, thermo-stability, and chemical resistance, all at a cost-effective price point.

**CAPP Reagent Bottles** come in a wide variety of sizes and volumes, with both wide and narrow mouth options making them perfect for storing liquids and solids alike. Made of the finest polyethylene and polypropylene, they offer excellent chemical resistance and thermal stability for both freezing and autoclaving; leaving you with little to worry for both shortand long-term storage applications.

Our narrow mouth reagent bottle sizes include 4, 8, 15, 30, 60, 125, 250 and 500 mL, while the wide mouth options are 30, 60, 125, 250 and 500 mL. These can be further configured with high density polyethylene (HDPE), low density polyethylene (LDPE) and polypropylene (PP) configurations. **CAPP Reagent Bottles** with LDPE and HDPE configurations can withstand temperatures from -100 to 120 degrees °C, making them freezer-friendly; while their PP counterparts are heat-stable from 0 to 135°C, perfect for autoclaving.

So, what are waiting for? Order your set of CAPP Reagent Bottles today! We know you will not regret it.



## Features:

- Versatile sizes and volumes with wide and narrow mouth options
- Chemically-resistant and leakproof for short- and long-term storage
- Made of high-quality LDPE/HDPE and autoclavable polypropylene
- Wide thermal stability range
- Made in a Class 100K clean room facility



## **Ordering Information**

Cat. No.	Description	Packaging	
PP Reagent Bottles, Narrow Mouth			
CAPP-PP-N-4	Bottle w/ cap, 4mL	6x12 pcs.	
CAPP-PP-N-8	Bottle w/ cap, 8mL	6x12 pcs.	
CAPP-PP-N-15	Bottle w/ cap, 15mL	6x12 pcs.	
CAPP-PP-N-30	Bottle w/ cap, 30mL	6x12 pcs.	
CAPP-PP-N-60	Bottle w/ cap, 60mL	6x12 pcs.	
CAPP-PP-N-125	Bottle w/ cap, 125mL	6x12 pcs.	
CAPP-PP-N-250	Bottle w/ cap, 250mL	6x12 pcs.	
CAPP-PP-N-500	Bottle w/ cap, 500mL	4x12 pcs.	
PP Reagent Bottles, Wide Mouth			
CAPP-PP-W-30	Bottle w/ cap, 30mL	6x12 pcs.	
CAPP-PP-W-60	Bottle w/ cap, 60mL	6x12 pcs.	
CAPP-PP-W-125	Bottle w/ cap, 125mL	6x12 pcs.	
CAPP-PP-W-250	Bottle w/ cap, 250mL	6x12 pcs.	
CAPP-PP-W-500	Bottle w/ cap, 500mL	4x12 pcs.	
HDPE Reagent Bottles, Narrow Mouth			
CAPP-HD-N-4	Bottle w/ cap, 4mL	6x12 pcs.	
CAPP-HD-N-8	Bottle w/ cap, 8mL	6x12 pcs.	
CAPP-HD-N-15	Bottle w/ cap, 15mL	6x12 pcs.	
CAPP-HD-N-30	Bottle w/ cap, 30mL	6x12 pcs.	
CAPP-HD-N-60	Bottle w/ cap, 60mL	6x12 pcs.	
CAPP-HD-N-125	Bottle w/ cap, 125mL	6x12 pcs.	
CAPP-HD-N-250	Bottle w/ cap, 250mL	6x12 pcs.	
CAPP-HD-N-500	Bottle w/ cap, 500mL	4x12 pcs.	
HDPE Reagent Bottles, Narrow Mouth, Amber Colour			
CAPP-HD-N-4-A	Bottle w/ cap, 4mL, Amber,	6x12 pcs.	
CAPP-HD-N-8-A	Bottle w/ cap 8mL, Amber,	6x12 pcs.	
CAPP-HD-N-15-A	Bottle w/ cap, 15mL, Amber	6x12 pcs.	
CAPP-HD-N-30-A	Bottle w/ cap, 30mL, Amber	6x12 pcs.	
CAPP-HD-N-60-A	Bottle w/ cap, 60mL, Amber	6x12 pcs.	
CAPP-HD-N-125-A	Bottle w/ cap 125mL, Amber,	6x12 pcs.	
CAPP-HD-N-250-A	Bottle w/ cap, 250mL, Amber	6x12 pcs.	
CAPP-HD-N-500-A	Bottle w/ cap, 500mL, Amber	4x12 pcs.	

Cat. No.	Description	Packaging	
HDPE Reagent Bottles, Wide Mouth			
CAPP-HD-W-30	Bottle w/ cap, 30mL	6x12 pcs.	
CAPP-HD-W-60	Bottle w/ cap, 60mL	6x12 pcs.	
CAPP-HD-W-125	Bottle w/ cap, 125mL	6x12 pcs.	
CAPP-HD-W-250	Bottle w/ cap, 250mL	6x12 pcs.	
CAPP-HD-W-500	Bottle w/ cap, 500mL	4x12 pcs.	
CAPP-HD-S-W-30	Bottle w/ cap, 30mL, sterile, Indiv. wrapped	500 pcs.	
CAPP-HD-S-W-60	Bottle w/ cap, 60mL, sterile, Indiv. wrapped	250 pcs.	
CAPP-HD-S-W-125	Bottle w/ cap, 125mL, sterile, Indiv. wrapped	125 pcs.	
CAPP-HD-S-W-250	Bottle w/ cap, 250mL, sterile, Indiv. wrapped	60 pcs.	
CAPP-HD-S-W-500	Bottle w/ cap, 500mL, sterile, Indiv. wrapped	30 pcs.	
HDPE Reagent Bottles, Wide Mouth, Amber Colour			
CAPP-HD-W-30-A	Bottle w/ cap, 30mL, Amber	6x12 pcs.	
CAPP-HD-W-60-A	Bottle w/ cap, 60mL, Amber	6x12 pcs.	
CAPP-HD-W-125-A	Bottle w/ cap, 125mL, Amber	6x12 pcs.	
CAPP-HD-W-250-A	Bottle w/ cap, 250mL, Amber	6x12 pcs.	
CAPP-HD-W-500-A	Bottle w/ cap, 500mL, Amber	4x12 pcs.	
LDPE Reagent Bottles, Narrow Mouth			
CAPP-LD-N-4	Bottle w/ cap, 4mL	6x12 pcs.	
CAPP-LD-N-8	Bottle w/ cap, 8mL	6x12 pcs.	
CAPP-LD-N-15	Bottle w/ cap, 15mL	6x12 pcs.	
CAPP-LD-N-30	Bottle w/ cap, 30mL	6x12 pcs.	
CAPP-LD-N-60	Bottle w/ cap, 60mL	6x12 pcs.	
CAPP-LD-N-125	Bottle w/ cap, 125mL	6x12 pcs.	
CAPP-LD-N-250	Bottle w/ cap, 250mL	6x12 pcs.	
CAPP-LD-N-500	Bottle w/ cap, 500mL	4x12 pcs.	
LDPE Reagent Bottles, Wide Mouth			
CAPP-LD-W-30	Bottle w/ cap, 30mL	6x12 pcs.	
CAPP-LD-W-60	Bottle w/ cap, 60mL	6x12 pcs.	
CAPP-LD-W-125	Bottle w/ cap, 125mL	6x12 pcs.	
CAPP-LD-W-250	Bottle w/ cap, 250mL	6x12 pcs.	
CAPP-LD-W-500	Bottle w/ cap, 500mL	4x12 pcs.	

