

## CSPP Harmony

SEROLOGICAL PIPETTES

CAPP Harmony Serological Pipettes are essential tools for accurate liquid measurement and transfer. With their precise graduation marks, colour-coded heads, and filter plugs, these pipettes ensure convenience, accuracy, and contamination-free experiments. Made from durable polystyrene, they are compatible with various pipettors and conform to USP CLASS VI standards. Trust CAPP Harmony Serological Pipettes to deliver reliable results in your laboratory or research facility.

## **Features:**

 Designed for accurate and rapid fluid dispensing, these polystyrene pipettes are sterile, pyrogen-free and cotton-plugged for a variety of liquid handling applications

• The sharp, black, highly legible permanent graduations prevent ambiguity and guarantee precise reading and accurate dispensing to +2%

 Pipettes are supplied individually wrapped in fibre-free colour-coded paper and plastic combination that is easy to open

• Each pipette comes with an effective filter barrier to prevent overflow of liquid into pipette

- Reverse graduations show the volume removed or remaining on all pipettes larger than 1mL
- Generous negative graduations provide additional pipetting volume
- The uniform mouth piece fits all major pipettors and is colour-coded for easy identification



## **Ordering Information**

Cat. No.	Description	Volume	Package	Colour
SP-1-C	CAPPHarmony Serological Pipette	1mL	individually wrapped, 500 pcs.	•
SP-2-C	CAPPHarmony Serological Pipette	2mL	individually wrapped, 500 pcs.	
SP-5-C	CAPPHarmony Serological Pipette	5mL	individually wrapped, 200 pcs.	
SP-10-C	CAPPHarmony Serological Pipette	10mL	individually wrapped, 200 pcs.	
SP-25-C	CAPPHarmony Serological Pipette	25mL	individually wrapped, 150 pcs.	
SP-50-C	CAPPHarmony Serological Pipette	50mL	individually wrapped, 100 pcs.	
SP-100-C	CAPPHarmony Serological Pipette	100mL	individually wrapped, 50 pcs.	•

