CAPP Microbiology
SINGLE CHANNEL PIPETTES

CAPP *Microbiology*. **Single Channel Pipettes** are ideal for microbiology tests in the food industry.

This pipette is supplied with three fixed volume knobs of 100μ L, 900μ L and 1000μ L and has a modified tip cone which is suited for use with a sterile extended length straw. This allows further reach and pipetting of viscous liquids such as milk, cream, yoghurt and cosmetics.

CAPP*Microbiology* **Single Channel Pipettes** are ideal instruments for protocols involving a dilution of 1:10.

Highlights:

- Enables extended reach
- Rotation-free volume selection
- Pre-set calibrated volume controllers
- Customisable volume selection
- Fully autoclavable



Ordering Information

CAPP Microbiology Pipettes (for use with sterile straws)

Cat. No.	Description	Inaccuracy %	Imprecision %
C1000-1BS	CAPPMicrobiology pipette, with 100μL, 900μL and 1000μL volume knobs	1.50/0.80	1.50/0.40
4170807	Sterile straws for CAPPMicrobiology pipette, 40x25 pcs.		



