

ELECTRONIC PIPETTES

CAPP



In addition to standard pipetting, you can choose between the most widely used pipetting modes, including reverse pipetting, multi dispensing, mixing and sequential dispensing.

Up to nine individual protocols can be stored, thanks to the extensive memory capacity.

Supplied with one long-life and easily exchangeable lithium battery, the **CAPP Transe Electronic Pipette** allows uninterrupted pipetting without any need to stop working in order to recharge.

The O-rings on the tip base provide an excellent fit, with a great variety of universal fitting tips.







Ordering Information

	Cat. No.	Description	Inaccuracy %	Imprecision %	Colour
	E20-8BA	CAPPTronic electronic pipette, 8-channel, 2-20μL	1.00	0.30	
	E200-8BA	CAPPTronic electronic pipette, 8-channel, 10-200μL	0.60	0.15	0
	E1200-8BA	CAPPTronic electronic pipette, 8-channel, 100-1200μL	0.50	0.15	

