

## REAGENT RESERVOIRS

**CAPP** *Origami* represents a new way of thinking when it comes to reagent reservoirs.

**CAPP***Origami* **Reagent Reservoirs** feature a sharp V-shaped bottom which helps prevent bubbles forming and ensures that you can retrieve the very last drop of liquid. This makes pipetting easier and ensures the most efficient use of the reagent.

The stiff construction of **CAPP** *Origami* results in a consistently firm reservoir. This easy-to-use design features pre-folded lines. To assemble your **CAPP** *Origami* **Reagent Reservoirs** simply flick it up by pushing the sides together, creating a strong, stable reservoir.

Made from durable polypropylene, these reservoirs are packed, shipped and stored as flat sheets, saving on both space and shipping costs.

Choose from 30mL reservoirs for 8/16 channel or 40mL for 12 channel pipettes.

Use your reagent most efficiently with CAPP Origami Reagent Reservoirs. Being environmentally and storage friendly, this product will increase space on your worktop and save you money.



## **Ordering Information**

Cat. No.	Description
CA40505	CAPPOrigami 30mL (8- and 16-channel pipettes), bag w/ 50 pcs.
CA40506	CAPPOrigami 30mL (8- and 16-channel pipettes), pre-sterile, 10 bags w/ 5 pcs. each
CA40510	CAPPOrigami 40mL (12-channel pipettes), bag w/ 50 pcs.
CA40511	CAPPOrigami 40mL (12-channel pipettes), pre-sterile, 10 bags w/ 5 pcs. each



## www.capp.dk • +49(0)3631/65242-0 • info@cappahn.com