

## **Highlights:**

- High quality metal shaft and piston assembly designed for greater tolerance to everyday laboratory conditions
- Ergonomic tip ejection lever with minimum thumb workload when offloading tips
- Stringent factory calibration protocol for each channel results in repeatable and dependable volume accuracy across all channels
- Comfortable in-hand feel narrows the possibility of repetitive strain injury





## **Ordering Information**

Cat. No.	Description	Inaccuracy %	Imprecision %	Colour
C02-16	CAPP multi pipettes, 16-channel, 0.2-2μL	4.00	2.00	
C10-16	CAPP multi pipettes, 16-channel, 0.5-10μL	2.00	1.00	$\circ$
C50-16	CAPP multi pipettes, 16-channel, 5-50μL	1.40	0.60	
C10-24	CAPP multi pipettes, 24-channel, 0.5-10µl	2.00	1.00	$\circ$
C50-24	CAPP multi pipettes, 24-channel, 5-50µl	1.40	0.60	
C10-32	CAPP multi pipettes, 32-channel, 0.5-10μl	4.00	2.00	$\bigcirc$
C10-48	CAPP multi pipettes, 48-channel, 0.5-10μL	4.00	2.00	$\circ$
C10-64	CAPP multi pipettes, 64-channel, 0.5-10μL	4.00	2.00	0

