CSPP Controller

MOTORISED PIPETTE CONTROLLER

Features:

• Light and fatigue free: weighs only 167g including the battery • Lithium battery with no memory effect and up to 12 hours of continuous use • Can be charged while in use • Powerful: 50mL pipette is filled in less than 10 seconds Express dispensing: depress the dispensing button beyond the built-in point of resistance Gravity delivery and blow-out dispensing mode selection Fully adjustable speed-setting • Exchangeable 0.2 μm PTFE disc filter acts as effective shut-down system if overfill occurs Exchangeable pipette holders are made of silicon • For use with glass or plastic pipettes ranging from 0.1mL to 100mL Vapour outlet for venting out corrosive liquids or fumes Fully autoclavable nose cone UV-resistant housing Complies with CE and RoHS

Ordering Information

Cat. No. Description PA-100 CAPPController Motorised Controller, 0.1-100mL

ddD



www.capp.dk • +49(0)3631/65242-0 • info@cappahn.com