Capp R10

MECHANICAL REPEATER

Features:

- Highly resistant towards a broad range of chemicals
- Easier selection and dispensing through a vision panel
- Reliable and maintenance-free due to its robust design
- Flexible: 12 tip sizes, 10 volume adjustments and 120 program steps
- Reliable ergonomics for your lab routine without sacrificing ultimate precision
- Handy, having a perfect ergonomic design. One-hand operation of volume selection, loading and dispensing of liquids
- Convenient: the easy-moving dosing button and the lightweight instrument makes long dosing series more comfortable
- High compatibility: **Capp** *R10* mechanical repeater is compatible both with **Capp** *Harmony* syringes and also with tips of many common brands
- The table below shows all dosing possibilities at a glance for choosing the most suitable tip for your application



Capp ₹10 Dispensing Volumes

					•	0				
Program No.	1	2	3	4	5	6	7	8	9	10
0.05mL	1 μL	1.5µL	2μL	2.5µL	3µL	3.5µL	4µL	4.5μL	5μL	5.5μL
0.1mL	2 μL	3µL	4µL	5μL	6μL	7μL	8µL	9μL	10μL	11μL
0.2mL	4 μL	6µL	8µL	10μL	12μL	14μL	16μL	18μL	20μL	22μL
0.5mL	10 μL	15μL	20μL	25μL	30μL	35μL	40μL	45μL	50μL	55μL
1mL	20 μL	30μL	40μL	50μL	60μL	70μL	80μL	90μL	100μL	110μL
1.25mL	25 μL	37.5μL	50μL	62.5μL	75μL	87.5μL	100μL	112.5μL	125μL	137.5μL
2.5mL	50 μL	75μL	100μL	125μL	150μL	175μL	200μL	225µL	250μL	275μL
5mL	100 μL	150μL	200μL	250μL	300μL	350μL	400μL	450μL	500μL	550μL
10mL	200µL	300µL	400μL	500μL	600μL	700μL	800µL	900µL	1000μL	1100μL
12.5mL	250µL	375μL	500μL	625μL	750µL	875μL	1000μL	1125µL	1250μL	1375μL
25mL	500 μL	750μL	1000μL	1250μL	1500μL	1750μL	2000μL	2250µL	2500μL	2750μL
50mL	1000 μL	1500μL	2000μL	2500μL	3000μL	3500μL	4000μL	4500μL	5000μL	5500μL
Steps	48	31	22	18	14	12	10	9	8	7

Ordering Information

Cat. No.	Description
PR-10	Capp R10 Mechanical Repeater

