

# CIPP Tempo PIPETTE CONTROLLER

The CAPP Tempo Pipette Controller has revolutionary features that separate it from other pipette controllers.

#### LED Display:

Display shows aspiration and dispense speeds
Displays remaining battery life and flashes when the power is low

- Displays remaining battery me and masters when the power is low
   Display turns-off after 20 seconds of non-use to conserve battery
  - Display turns-off after 20 seconds of non-use to conserve battery

### Adjustable Speeds:

Powerful motor fills 25mL in 5 seconds (at speed level 6)
Aspiration and dispense modes have separate speed controls
Six different speed settings for both aspiration and dispense modes

### Electronic Brake System (EBS):

- Apply EBS to slow down at the last second to dial in the most accurate volume
- Press EBS a second time to return to maximum speed
- EBS maximises speed and precision

### **Power Supply:**

• Lithium-ion battery, full charge provides eight hours of continuous use

the the the

• Battery can be charged while pipette is in use

## **Ordering Information**

Cat. No.	Description
Т100-В	CAPPTempo Blue, 0.1-100mL
T100-G	CAPPTempo Green, 0.1-100mL
T100-R	CAPPTempo Red, 0.1-100mL
T100-V	CAPPTempo Violet, 0.1-100mL
T100-CS	CAPPTempo Charging Stand



111 AUDIT