EXPell Cryoboxes

CARDBOARD CRYOBOXES

These **EXPell** Cryolores from **CAPP** are made out of high-quality cardboard and are available with two different grid sizes: 9x9 (81 cells) and 10x10 (100 cells). For easy sample detection our boxes feature a write-on surface which will make identifying and retrieving your samples a lot easier. **CAPP**'s cardboard **EXPell** Cryolores are chemically-resistant to alcohols and mild organic solvents, but you cannot autoclave them.

If you are searching for an environmentally-friendly option for your cryo storage these are a great option as they are fully disposable after use.

Features:

- Cardboard is constructed with a water-repellent finish for durability
- Temperature-resistant, recyclable cardboard
- Removable grid allowing more flexibility for your needs
- Supplied in 9x9 and 10x10 grid options

With our cryoboxes you don't have to sacrifice quality to save the environment.

Both are important!

Ordering Information

Cat. No.	Description
5042030C	Cardboard Cryobox B9, 133 x 133 x 50mm, 9x9 cell grid, case with 36 pcs.
5042020C	Cardboard Cryobox B10, 133 x 133 x 50mm including 10x10 cell grid, case with 36 pcs.

For more information about other size options, please contact info@cappahn.com

