

PCR MICROTUBES, STRIPS & PLATES OF 0.2ML (REGULAR PROFILE)

Experience reliable PCR results with CAPP EXPell *PCR* Microtubes, Strips and Plates.

Made from high-quality polypropylene, these consumables are free of contaminants, ensuring accurate and consistent performance. The alphanumeric coding and raised

well rim of the PCR Plates facilitate easy sample tracking and secure sealing, while frosted and white options enhance signal-to-noise ratios.

With multiple skirt options and compatibility with standard and Real-Time PCR routines, **CAPP EXPell** *PCR*, consumables are the trusted choice for your PCR needs. You can rely on **CAPP** for purity, reliability, and precise PCR results.

CAPP also provides a variety of sealing adhesives for optimal closure of the plates and other accessories.

Please refer to the ordering information.

Ordering Information

Cat. No.	Description
5100200C	EXPellPlus PCR tubes 0.2mL, bag, 20x1000 pcs.
5100201C	EXPell PCR tubes 0.2mL, bag, 40x500 pcs.
5100205C	EXPell PCR tubes 0.2mL, bag, pre-sterile 40x500 pcs.
5100210C	EXPell PCR 8-strip and separate flat caps 0.2mL, 5x10x12 pcs.
5100212C	EXPell PCR 8-strip w/ single attached flat caps 0.2mL, 5x10x12 pcs.
5100217C	EXPell PCR 8-strip and separate optical caps 0.2mL, 5x10x12 pcs.
5100218C	EXPell PCR 8-strip, white, and separate optical caps 0.2mL, 5x10x12 pcs.
5100419C	Expell PCR 96-well plate, frosted, non-skirted, 0.2mL, 5x10 pcs.
5100420C	Expell PCR 96-well plate, frosted, sub semi-skirted, 0.2mL, 5x10 pcs.
5100421C	Expell PCR 96-well plate, frosted, semi-skirted, 0.2mL, 5x10 pcs.

