

PIPETTE STARTER KITS

CAPP offers a range of convenient Pipette Starter Kits, which contain three single channel pipettes, three racks of corresponding tips, a pipette stand and a handy laboratory timer.

Made of durable materials and quality control tested, the pipettes and tips in our Pipette Starter Kits offer accuracy and precision like no other.

Choosing does not have to be a daunting task with Pipette Starter Kits from CAPP! Our versatile offerings are sure to meet the needs of your laboratory, no matter the intended application. Our kits offer quality, durability and industry-standard accuracy while providing a variety of options and models.

There are three pipette model kits to choose from:







Features:

- Included in each Pipette Starter Kit are our Expell and ExpellPlus racked tips which are made of low retention virgin polypropylene and are certified RNase, DNase and Pyrogen free
- All pipettes and tips are fully autoclavable and meet industry standards for precision and accuracy
- Each Starter Kit contains a pipette stand and laboratory timer for facilitating your work
- Our compact kits come conveniently packaged in space-saving, recyclable cardboard boxes for easy storage
- Users can choose from the following options:
 - 10 mL, 100 mL and 1000 mL volume kits
 - 20 mL, 200 mL and 1000 mL volume kits
 - Any three sizes of pipette from their model of choice (CAPPBravo and ecopipette only)



Designed for academic, diagnostic and start-up laboratories with limited budgets, the CAPPSolo Starter Kit is a great bargain for budget-conscious buyers who still want a precise and accurate product for their pipetting applications.







For the user seeking a high quality, ergonomic experience, these lightweight CAPPBravo pipettes make this Pipette Starter Kit an excellent pick. Like all CAPP pipettes, the CAPPBravo line is compatible with CAPP EXPellPlus pipette tips and most other major tip brands.





The ecopipette Starter Kit is a great option for the eco-conscious buyer who is looking for a dependable, accurate pipette set.

Packaged using easily recyclable eco-friendly materials, the ecopipette is constructed of some of the most renewable resources available, while still maintaining quality and durability.



Think SUSTAINABLY!

REDUCE the plastic, metal and rubber waste at your laboratory with the ecotrade programme.

Contact <u>info@cappahn.com</u> for more information.

03

ENVIRONMENTALLY-FRIENDLY RECYCLING

As part of the Green Dot scheme (Der Grüne Punkt), AHN Biotechnologie GmbH will ensure the environmentally friendly recycling of all components



Recycle

01

PIPETTE RETURN

Send your old pipettes, regardless of brand and condition, to your local distributor who will forward them to AHN Biotechnologie GmbH in Nordhausen, Germany

02

DISASSEMBLY AND SORTING

Our dedicated team in Nordhausen will disassemble all the components and sort them into plastic, rubber, metal, etc., in accordance with the German national recycling programme

Ordering Information

ecotrade

CAPP Solo Starter Kit

SK-03-XX

Cat. No.	Description
SK-03-S01	Three CAPPSolo pipettes S10-1 /S100-1 /S1000-1, one C-05 stand, three racks of tips and one CAPPtimer
SK-03-S02	Three CAPPSolo pipettes S20-1 /S200-1 /S1000-1, one C-05 stand, three racks of tips and one CAPPtimer
CAPP Bravo Starter Kit	
Cat. No.	Description
SK-03-B01	Three CAPPBravo pipettes B10-1 /B100-1 /B1000-1, one C-05 stand, three racks of tips and one CAPPtimer
SK-03-B02	Three CAPPBravo pipettes B20-1 /B200-1 /B1000-1, one C-05 stand, three racks of tips and one CAPPtimer
SK-03-BB	Three CAPPBravo pipettes of choice, one C-05 stand, three racks of tips and one CAPPtimer
ecopipette Starter Kit	
Cat. No.	Description
SK-03-E01	Three ecopipettes C10-1 /C100-1 /C1000-1, one C-05 stand, three racks of tips and one CAPPtimer
SK-03-E02	Three ecopipettes C20-1 /C200-1 /C1000-1, one C-05 stand, three racks of tips and one CAPPtimer

Three ecopipettes of choice, one C-01 stand, three racks of tips and one CAPPtimer