

CAPP *Harmony*

COMPATIBILITY CHART

For syringes' compatibility with other pipette brands, please contact info@capp.dk



CAPP *Harmony* Basic Retips

RS-01-C	RS-03-C	RS-06-C
RS-01S-C	RS-03S-C	RS-06S-C
0.05mL 100 pcs	1.25mL 100 pcs	12.5mL 100 pcs
Fully compatible	Fully compatible	Fully compatible
Fully compatible	Fully compatible	Fully compatible
Fully compatible	Fully compatible	Fully compatible
Not compatible	Not compatible	Not compatible
Not compatible	Not compatible	Not compatible
Fully compatible	Fully compatible	Fully compatible
Not compatible	Not compatible	Not compatible
Fully compatible	Fully compatible	Fully compatible
Fully compatible	Fully compatible	Fully compatible
Fully compatible	Fully compatible	Fully compatible
Not compatible	Fully compatible	Fully compatible
Fully compatible	Fully compatible	Fully compatible
Fully compatible	Fully compatible	Fully compatible

CAPP *Harmony* Macro Retips

RS-07-C	RS-08-C
RS-07S-C	RS-08S-C
25mL 25 pcs	50mL 25 pcs
Fully compatible	Fully compatible
Fully compatible	Fully compatible
Fully compatible	Fully compatible
Fully compatible	Fully compatible
Fully compatible	Fully compatible
Fully compatible	Fully compatible
Fully compatible	Fully compatible
Fully compatible	Fully compatible
Fully compatible	Fully compatible
Fully compatible	Fully compatible
Not compatible	Not compatible
Fully compatible	Fully compatible
Fully compatible	Fully compatible

CAPP Rhythm Repeater
CAPP R10 Repeater
Biohit Mechanical Stepper
eLINE Lite Dispenser
eLINE Pro Dispenser
Fisherbrand
Eppendorf Multipipette Plus
HTL MINILAB 201
Brand HandyStep
Brand HandyStep Electronic
Gilson Distriman / Repetman
Raining Auto Rep S
Raining Auto Rep E

© 2015 CAPP ApS. All rights reserved, CAPP Bravo, ecopipette and CAPPAero are trademarks of CAPP ApS. eline, mline and Proline Plus are trademarks of Sartorius Biohit Liquid Handling Oy. Transferpette is a trademark of BrandTech Scientific Inc. Eppendorf Reference and Eppendorf Research Plus are registered trademarks of Eppendorf AG. Finnipipette is a trademark of Thermo Fisher Scientific Inc. Pipetman is a trademark of Gilson, Inc. Rainin Classic is a trademark of Mettler-Toledo International, Inc. WJR brand is a trademark of WJR International, Inc. Labmate is a trademark of HTLPZ HTL SA. Biopette is a trademark of Labnet International, Inc.

CAPP *Harmony*

COMPATIBILITY CHART

For syringes' compatibility with other pipette brands, please contact info@capp.dk



CAPP *Harmony* Universal Retips

RS-11-C RS-11S-C 0.1mL 100 pcs	RS-12-C RS-12S-C 0.2mL 100 pcs	RS-02-C RS-02S-C 0.5mL 100 pcs	RS-13-C RS-13S-C 1mL 100 pcs	RS-04-C RS-04S-C 2.5mL 100 pcs	RS-05-C RS-05S-C 5mL 100 pcs	RS-14-C RS-14S-C 10mL 100 pcs
Not compatible	Not compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible	Not compatible
Fully compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible
Not compatible	Not compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible	Not compatible
Fully compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible
Fully compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible
Not compatible	Not compatible	Fully compatible	Not compatible	Fully compatible	Fully compatible	Not compatible
Fully compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible
Not compatible	Not compatible	Fully compatible	Not compatible	Fully compatible	Fully compatible	Not compatible
Fully compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible
Fully compatible	Fully compatible	Fully compatible	Not compatible	Fully compatible	Fully compatible	Not compatible
Not compatible	Not compatible	Not compatible	Not compatible	Not compatible	Not compatible	Not compatible
Fully compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible
Fully compatible	Fully compatible	Fully compatible	Not compatible	Fully compatible	Fully compatible	Not compatible

© 2015 CAPP ApS. All rights reserved, CAPP Bravo, ecopipette and CAPPAero are trademarks of CAPP ApS. eline, mline and Proline Plus are trademarks of Sartorius Biohit Liquid Handling Oy. Transferpette is a trademark of BrandTech Scientific Inc. Eppendorf Reference and Eppendorf Research Plus are registered trademarks of Eppendorf AG. Finnpiquette is a trademark of Thermo Fisher Scientific Inc. Pipetman is a trademark of Gilson, Inc. Rainin Classic is a trademark of Mettler-Toledo International, Inc. WJR brand is a trademark of WJR International, Inc. Labmate is a trademark of HTLPZ HTL SA. Biopette is a trademark of Labnet International, Inc.

