

Capp Harmony

COMPATIBILITY CHART

For syringes compatibility with other pipette brands, please contact info@capp.dk



**RS-01-C
RS-01S-C**
Capp Harmony
syringes 0,05 ml,
100 pcs.



**RS-11-C
RS-11S-C**
Capp Harmony
syringes 0,1 ml, 100
pcs.



**RS-12-C
RS-12S-C**
Capp Harmony
syringes 0,2 ml,
100 pcs.



**RS-02-C
RS-02S-C**
Capp Harmony
syringes 0,5 ml,
100 pcs.



**RS-13-C
RS-13S-C**
Capp Harmony
syringes 1,0 ml,
100 pcs.




**RS-03-C
RS-03S-C**
Capp Harmony
syringes 1,25 ml,
100 pcs.

Capp Rhythm repeater	Fully compatible	Not compatible	Not compatible	Fully compatible	Not compatible	Fully compatible
Capp R10 repeater	Fully compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible
Biohit Mechanical Stepper	Fully compatible	Not compatible	Not compatible	Fully compatible	Fully compatible	Fully compatible
eLINE Lite Dispenser	Not compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible	Not compatible
eLINE Pro Dispenser	Not compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible	Not compatible
Fisherbrand	Fully compatible	Not compatible	Not compatible	Fully compatible	Not compatible	Fully compatible
HTL MINILAB 201	Fully compatible	Not compatible	Not compatible	Fully compatible	Not compatible	Fully compatible
Brand HandyStep	Fully compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible
Brand HandyStep electronic	Fully compatible	Fully compatible	Fully compatible	Fully compatible	Not compatible	Fully compatible
Gilson Distriman/ Repetman	Not compatible	Not compatible	Not compatible	Not compatible	Not compatible	Fully compatible
Rainin AutoRep S	Fully compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible
Rainin AutoRep E	Fully compatible	Fully compatible	Fully compatible	Fully compatible	Not compatible	Fully compatible

TEL. (+45) 6613 6140
 FAX. (+45) 6613 2770

MAIL: INFO@CAPP.DK
 WEB: WWW.CAPP.DK



Fully compatible	Fully compatible	Not compatible	Fully compatible	Fully compatible	Fully compatible
Fully compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible
Fully compatible	Fully compatible	Not compatible	Fully compatible	Fully compatible	Fully compatible
Fully compatible	Fully compatible	Fully compatible	Not compatible	Fully compatible	Fully compatible
Fully compatible	Fully compatible	Fully compatible	Not compatible	Fully compatible	Fully compatible
Fully compatible	Fully compatible	Not compatible	Fully compatible	Fully compatible	Fully compatible
Fully compatible	Fully compatible	Not compatible	Fully compatible	Fully compatible	Fully compatible
Fully compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible
Fully compatible	Fully compatible	Not compatible	Fully compatible	Fully compatible	Fully compatible
Not compatible	Not compatible	Not compatible	Fully compatible	Not compatible	Not compatible
Fully compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible	Fully compatible
Fully compatible	Fully compatible	Not compatible	Fully compatible	Fully compatible	Fully compatible