

EXPELLPLUS

INDIVIDUALLY WRAPPED TIPS

EXPELLPLUS individually wrapped tips are ideal for sensitive applications requiring the strictest aseptic conditions.

Their unique individual packing ensures sterility for every single pipetting cycle, maximizing the exceptional performance of EXPELLPLUS sterile tips. Assuring the greatest purity standards, EXPELLPLUS individually wrapped tips are free from Endotoxin, DNA, DNase, RNase and PCR inhibitor.

Available in four volume sizes: 10µL standard and extended length, 200µL and 1000-1250µL, EXPELLPLUS individually wrapped tips are universally compatible with the majority of pipette brands available on the market. Thanks to the unique composition of plastic polymer, the inner surface of EXPELLPLUS individually wrapped tips is more hydrophobic, which significantly reduces sample adhesion. These low retention properties minimize sample binding inside the tip, making EXPELLPLUS tips ideal for applications requiring high accuracy and reproducibility.

EXPELLPLUS individually wrapped tips are sterile. Batch number and expiration date are printed on each individual wrap.



Ordering Information

Cat. No. Description

5030012C	ExpellPlus 10µL, pre-sterile, clear, bag, 5x400 pcs. individually wrapped tips
5030042C	ExpellPlus 10µL XL, pre-sterile, long, clear, bag, 5x400 pcs. individually wrapped tips
5030072C	ExpellPlus 200µL, pre-sterile, clear, bag, 5x400 pcs. individually wrapped tips
5130132C	Expell 1000µL (1250µL), pre-sterile, clear, bag, 5x200 pcs. individually wrapped tips

www.capp.dk – +49 (0) 3631 65242-0 – info@capp.dk