

Expell

BLACK MICROCENTRIFUGE TUBES

With a special focus on processing light-sensitive samples, Capp has developed **Expell** black microcentrifuge tubes. The black tubes are available in one size (1.5mL) and come bulk packed in bags of 500pcs.

The black Polypropylene compound creates a completely opaque barrier, preventing light from causing any reactions within the liquid inside the tube, making the tubes ideal for processing light sensitive samples. **Expell** black microcentrifuge tubes are a great alternative to the commonly used method of wrapping the tube in aluminum foil in order to protect the samples from light impact.

Expell black microcentrifuge tubes are suitable for applications requiring photo-sensitive, fluorescent or ethidium bromide probes. Unlike amber tinted tubes, **Expell** black microcentrifuge tubes are completely opaque and protect even the most light-sensitive samples.



Ordering Information

Expell Black Microcentrifuge Tubes

Cat. No.	Description
5101550C	Expell Microcentrifuge tubes 1.5 mL, bulk, black, 20x500 pcs