

# CAPP Solo

MULTICHANNEL PIPETTES

**CAPP Solo Multichannel Pipettes** offer industry-standard pipetting accuracy and precision in a very affordable package. Manufactured with budget-constrained customers in mind, these pipettes offer quick, reproducible pipetting without compromising on quality.

**CAPP Solo Multichannel Pipettes** are supplied in 8-channel and 12-channel variants and are guaranteed to last in academic, diagnostic and other comparable laboratories who require workhorses that can withstand extended usage cycles.

#### Features:

- Exchangeable colour inserts for easy identification of pipettes by user or volume capacities
- Remarkably lightweight pipette chassis for its class results in more comfortable pipette handling
- Counter-mechanism prevents accidental rotation while pipetting
- Accurate liquid dispensing across all channels for repeatable results
- Durable build quality ensures chemical resistance, bump protection and thermal stability
- Long-term calibration retention, even after sterilising in an autoclave
- Each tip cone is universally compatible with tips from the industry's most common brands
- Supplied with a standard two-year warranty





## Ordering Information

### CAPP Solo Multichannel Pipettes

Cat. No.	Description	Inaccuracy %	Imprecision %
S10-8	CAPPSolo 8-channel pipette, 0.5 – 10 $\mu$ L	5.00 / 2.00	3.60 / 1.00
S100-8	CAPPSolo 8-channel pipette, 10 – 100 $\mu$ L	2.50 / 1.30	1.20 / 0.50
S300-8	CAPPSolo 8-channel pipette, 30 – 300 $\mu$ L	1.66 / 1.40	0.60 / 0.30
S10-12	CAPPSolo 12-channel pipette, 0.5 – 10 $\mu$ L	5.00 / 2.00	3.60 / 1.00
S100-12	CAPPSolo 12-channel pipette, 10 – 100 $\mu$ L	2.50 / 1.30	1.20 / 0.50
S300-12	CAPPSolo 12-channel pipette, 30 – 300 $\mu$ L	1.66 / 1.40	0.60 / 0.30

