

# EX Pell Cryo Boxes

## CARDBOARD CRYOBOXES

These **EX Pell Cryoboxes** from **CAPP** are made out of high quality cardboard and are available with two different grid sizes; 9x9 (81 cells) and 10x10 (100 cells). For easy sample detection our boxes feature a write-on surface, which will make identifying and retrieving your samples a lot easier. **CAPP's** cardboard **EX Pell Cryoboxes** are chemically-resistant to alcohols and mild organic solvents, but you cannot autoclave them.

If you are searching for an environmental friendly option for your cryo storage these are a great option as they are fully disposable after use.

### Features:

- Cardboard is constructed with a water-repellent finish for durability
- Temperature resistant, recyclable cardboard
- Removable grid, allowing more flexibility for your needs
- Supplied as 9x9 and 10x10 grid options



## Ordering Information

Cat. No.	Description
5042030C	Cardboard Cryobox B9, 133x133x50 mm, 9x9 cell grid, case with 36 pcs.
5042020C	Cardboard Cryobox B10, 133x133x50 mm including 10x10 cell grid, case with 36 pcs.

For more information about other size options, please contact [info@cappahn.com](mailto:info@cappahn.com)

