



Expell PCR Plates Compatibility Chart
Thermocycler Compatibility for 5100214C
 ref. PIP_DS_5100214C_22032019

Abbott	m2000rt	96 x 0.2ml		Compatible
Agilent	SureCycler 8800 thermal cycler	96 x 0.2ml		Compatible
Applied Biosystems	310	96 x 0.2ml		Compatible
Applied Biosystems	3100	96 x 0.2ml		Compatible
Applied Biosystems	3130	96 x 0.2ml		Compatible
Applied Biosystems	3500 Genetic Analyzer	8-capillary		Compatible
Applied Biosystems	3700	96 x 0.2ml		Compatible
Applied Biosystems	3730	96 x 0.2ml		Compatible
Applied Biosystems	7300 Real-Time PCR System	96 x 0.2ml		Compatible
Applied Biosystems	7500 Real-Time PCR System	96 x 0.2ml Standard		Compatible
Applied Biosystems	7900HT Fast Real-Time PCR System	96 x 0.2ml Standard		Compatible
Applied Biosystems	ABI PRISM 7000 SDS	96 x 0.2ml		Compatible
Applied Biosystems	ABI PRISM 7700 SDS	96 x 0.2ml		Compatible
Applied Biosystems	GeneAmp PCR System 2700	96 x 0.2ml		Compatible
Applied Biosystems	GeneAmp PCR System 2720	96 x 0.2ml		Compatible
Applied Biosystems	GeneAmp PCR System 5700	96 x 0.2ml		Compatible
Applied Biosystems	GeneAmp PCR System 9600	96 x 0.2ml		Compatible
Applied Biosystems	GeneAmp PCR System 9700	96 x 0.2ml		Compatible
Applied Biosystems	GeneAmp PCR System 9700	Dual 96-Well 96 x 0.2ml		Compatible
Applied Biosystems	GeneAmp PCR System 9700 Aluminum	96 x 0.2ml		Compatible
Applied Biosystems	GeneAmp PCR System 9700 Gold	96 x 0.2ml		Compatible
Applied Biosystems	MiniAmp Plus Thermal Cycler	96 x 0.2ml		Compatible
Applied Biosystems	MiniAmp Thermal Cycler	96 x 0.2ml		Compatible
Applied Biosystems	ProFlex™ PCR System	96 x 0.2ml		Compatible
Applied Biosystems	QuantStudio™ 12k Flex	96 x 0.2ml		Compatible
Applied Biosystems	QuantStudio™ 3 Real-Time PCR System	96 x 0.2ml		Compatible
Applied Biosystems	QuantStudio™ 5 Real-Time PCR System	96 x 0.2ml		Compatible
Applied Biosystems	QuantStudio™ 6 Flex	96 x 0.2ml		Compatible
Applied Biosystems	QuantStudio™ 7 Flex	96 x 0.2ml		Compatible
Applied Biosystems	SimpliAmp™	96 x 0.2ml		Compatible
Applied Biosystems	Solid3	96 x 0.2ml		Compatible
Applied Biosystems	Veriti™ 96-Well Thermal Cycler	96 x 0.2ml		Compatible
Applied Biosystems	Veriti™ Dx 96-Well Thermal Cycler	96 x 0.2ml		Compatible
Applied Biosystems	ViiA™ 7 Real-Time PCR System	96 x 0.2ml		Compatible
Bioneer	MyGenie 96 Gradient Thermal Block	96 x 0.2ml		Compatible
DuPont Qualicon	BAX System Q7	96 x 0.2ml		Compatible
Siemens	VERSANT kPCR Molecular System AD	96 x 0.2ml		Compatible
Stratagene	Mx3000P™	96 x 0.2ml		Compatible
Stratagene	Mx3005P™	96 x 0.2ml		Compatible