



Expell PCR Plates Compatibility Chart
Thermocycler Compatibility for 5100314C
ref. PIP_DS_5100314C_22032019

nalytik-Jena	Alpha SC	96 x 0.2ml	✓	Compatible
Analytik-Jena	FlexCycler	Monoblock 96 x 0.2ml	✓	Compatible
Analytik-Jena - Biometra	TOne 96	96 x 0.2ml	✓	Compatible
Analytik-Jena - Biometra	TOne 96G	96 x 0.2ml	✓	Compatible
Apollo	ATC201	96-48 96 x 0.2ml	✓	Compatible
Apollo	ATC401	96 x 0.2ml	✓	Compatible
Apollo	ATC701	96 x 0.2ml	✓	Compatible
Apollo	ATC801	96 x 0.2ml	✓	Compatible
Axygen®	MaxyGene™ II	96 x 0.2ml	✓	Compatible
Barnstead Thermolyne	Amplitrion I	96 x 0.2ml	✓	Compatible
Barnstead Thermolyne	Amplitrion II	96 x 0.2ml	✓	Compatible
Bioer	GenePro B-96G	96 x 0.2ml	✓	Compatible
Bioer	Life Express TC-96/T/H(a)	96 x 0.2ml	✓	Compatible
Bioer	LifePro TC-96/G/H(b)A	96 x 0.2ml	✓	Compatible
Bio-Rad Laboratories	C1000	96-Deep Well	✓	Compatible
Bio-Rad Laboratories	C1000	96-Well Fast	✓	Compatible
Bio-Rad Laboratories	C1000 Touch	96-Deep Well	✓	Compatible
Bio-Rad Laboratories	C1000 Touch	96-Well Fast	✓	Compatible
Bio-Rad Laboratories	CFX96 Connect Real-Time PCR Detection System	96 x 0.2ml	✓	Compatible
Bio-Rad Laboratories	CFX96 Connect Real-Time PCR Detection System	96-Deep Well	✓	Compatible
Bio-Rad Laboratories	CFX96 Real-Time PCR Detection System	96 x 0.2ml	✓	Compatible
Bio-Rad Laboratories	CFX96 Real-Time PCR Detection System	96-Deep Well	✓	Compatible
Bio-Rad Laboratories	DNA Engine	96v Alpha unit 96 x 0.2ml	✓	Compatible
Bio-Rad Laboratories	DNA Engine	Chromo4 96 x 0.2ml	✓	Compatible
Bio-Rad Laboratories	DNA Engine Dyad	96v Alpha unit 96 x 0.2ml	✓	Compatible
Bio-Rad Laboratories	DNA Engine Dyad	Chromo4 96 x 0.2ml	✓	Compatible
Bio-Rad Laboratories	DNA Engine Dyad Disciple	96v Alpha unit 96 x 0.2ml	✓	Compatible
Bio-Rad Laboratories	DNA Engine Dyad Disciple	Chromo4 96 x 0.2ml	✓	Compatible
Bio-Rad Laboratories	DNA Engine Tetrad 2	96v Alpha unit 96 x 0.2ml	✓	Compatible
Bio-Rad Laboratories	DNA Engine Tetrad 2	Chromo4 96 x 0.2ml	✓	Compatible
Bio-Rad Laboratories	Opticon 2	96 x 0.2ml	✓	Compatible
Bio-Rad Laboratories	S1000	96-Deep Well	✓	Compatible
Bio-Rad Laboratories	S1000	96-Well Fast	✓	Compatible
Bio-Rad Laboratories	T100	96 x 0.2ml	✓	Compatible
Blue-Ray Biotech	TurboCycler TCST-9610	96 x 0.2ml	✓	Compatible
Blue-Ray Biotech	TurboCycler TCST-9620	96 x 0.2ml	✓	Compatible
Bulldog Bio	GenePro™ Thermal Cycler	96 x 0.2ml	✓	Compatible
Bulldog Bio	GeneTouch™ Thermal Cycler	96 x 0.2ml	✓	Compatible
Bulldog Bio	LifePro™ Thermal Cycler	96 x 0.2ml	✓	Compatible
Bulldog Bio	LifeTouch™ Thermal Cycler	96 x 0.2ml	✓	Compatible
eENZYME	GenePro Thermal Cycler	96 x 0.2ml	✓	Compatible
eENZYME	LifeExpress B-96D	96 x 0.2ml	✓	Compatible
eENZYME	LifeExpress TC-96	96 x 0.2ml	✓	Compatible
eENZYME	MyGene MG-96 Gradient Thermal Cycler	96-77 96 x 0.2ml	✓	Compatible
Eppendorf	Mastercycler	96-77 96 x 0.2ml	✓	Compatible
Eppendorf	Mastercycler ep gradient	96 x 0.2ml	✓	Compatible
Eppendorf	Mastercycler ep gradient Silver	96 x 0.2ml	✓	Compatible
Eppendorf	Mastercycler ep realplex	96 x 0.2ml	✓	Compatible
Eppendorf	Thermomixer C	96 x 0.2ml	✓	Compatible
Finnzymes Instruments	TCA0096 Arktik	96 x 0.2ml	✓	Compatible
Gene Technologies Ltd.	G-Storm GS1	96 x 0.2ml	✓	Compatible
Gene Technologies Ltd.	G-Storm GS1	96-48 96 x 0.2ml	✓	Compatible
Gene Technologies Ltd.	G-Storm GS2	96 x 0.2ml	✓	Compatible
Gene Technologies Ltd.	G-Storm GS4	96 x 0.2ml	✓	Compatible
Gene Technologies Ltd.	G-Storm GS4	96-48 96 x 0.2ml	✓	Compatible
Gene Technologies Ltd.	G-Storm GSX	96 x 0.2ml	✓	Compatible
Gene Technologies Ltd.	G-Storm GSX	96-48 96 x 0.2ml	✓	Compatible
Hain-LIFESCIENCE	GTQ-Cycler 96	96 x 0.2ml	✓	Compatible
Kyratec	SuperCycler	96 x 0.2ml	✓	Compatible
LongGene Scientific Instruments Co.,Ltd.	A100	96 x 0.2ml	✓	Compatible
LongGene Scientific Instruments Co.,Ltd.	A200	96-77 77 x 0.5ml	✓	Compatible
LongGene Scientific Instruments Co.,Ltd.	A200	96-77 96 x 0.2ml	✓	Compatible
MJ Research Inc.	Opticon 2 System	96 x 0.2ml	✓	fits cycler in whole or in part
MJ Research Inc.	Opticon System	96 x 0.2ml	✓	fits cycler in whole or in part
MJ Research Inc.	PTC-100	96 x 0.2ml	✓	Compatible
MJ Research Inc.	PTC-100	96 x 0.2ml AgV	✓	Compatible
MJ Research Inc.	PTC-200 DNA Engine	96v Alpha unit 96 x 0.2ml	✓	Compatible
MJ Research Inc.	PTC-200 DNA Engine	96v Moto Alpha Unit 96 x 0.2ml	✓	Compatible
MJ Research Inc.	PTC-200 DNA Engine	Chromo4 96 x 0.2ml	✓	Compatible
MJ Research Inc.	PTC-220 DNA Engine Dyad	96v Alpha unit 96 x 0.2ml	✓	Compatible
MJ Research Inc.	PTC-220 DNA Engine Dyad	96v Moto Alpha Unit 96 x 0.2ml	✓	Compatible
MJ Research Inc.	PTC-220 DNA Engine Dyad	Chromo4 96 x 0.2ml	✓	Compatible
MJ Research Inc.	PTC-221 DNA Engine Disciple	96v Alpha unit 96 x 0.2ml	✓	Compatible
MJ Research Inc.	PTC-221 DNA Engine Disciple	96v Moto Alpha Unit 96 x 0.2ml	✓	Compatible
MJ Research Inc.	PTC-221 DNA Engine Disciple	Chromo4 96 x 0.2ml	✓	Compatible
MJ Research Inc.	PTC-225 DNA Engine Tetrad	96v Alpha unit 96 x 0.2ml	✓	Compatible
MJ Research Inc.	PTC-225 DNA Engine Tetrad	96v Moto Alpha Unit 96 x 0.2ml	✓	Compatible

MJ Research Inc.	PTC-225 DNA Engine Tetrad	Chromo4 96 x 0.2ml	✓	Compatible
MJ Research Inc.	PTC-240 DNA Engine Tetrad 2	96v Alpha unit 96 x 0.2ml	✓	Compatible
MJ Research Inc.	PTC-240 DNA Engine Tetrad 2	96v Moto Alpha Unit 96 x 0.2ml	✓	Compatible
MJ Research Inc.	PTC-240 DNA Engine Tetrad 2	Chromo4 96 x 0.2ml	✓	Compatible
PCRMax	AlphaCycler AC-1	96 x 0.2ml	✓	Compatible
PCRMax	AlphaCycler AC-4	96 x 0.2ml	✓	Compatible
Peqlab	Primus 96 Advanced	96-77 96 x 0.2ml	✓	Compatible
Roche	cobas z 480	96 x 0.2ml	✓	Compatible
Roche	LightCycler @ 480 Real-Time PCR System	96 x 0.2ml	✓	Compatible
Roche	LightCycler @ 96 Real-Time PCR System	96 x 0.2ml	✓	Compatible
SensoQuest	Labcyler	Thermoblock 96	✓	Compatible
Techne	Flexigene	96 x 0.2ml	✓	Compatible
Techne	Gene-E	96 x 0.2ml	✓	Compatible
Techne	Genius	96 x 0.2ml	✓	Compatible
Techne	Prime thermal cycler	96 x 0.2ml	✓	Compatible
Techne	TC-4000	96 x 0.2ml	✓	Compatible
Techne	TC-412	96 x 0.2ml	✓	Compatible
Techne	TC-5000	96 x 0.2ml	✓	Compatible
Techne	TC-512	96 x 0.2ml	✓	Compatible
Techne	TC-Plus	96 x 0.2ml	✓	Compatible
Techne	Touchgene Gradient	96 x 0.2ml	✓	Compatible
ThermoFisher Scientific	Arktik Thermal Cycler	96 x 0.2ml	✓	Compatible
VWR	DuoCycler	96 x 0.2ml	✓	Compatible
VWR	DuoCycler fast	96 x 0.2ml	✓	Compatible
VWR	QuattroCycler	96 x 0.2ml	✓	Compatible
VWR	QuattroCycler fast	96 x 0.2ml	✓	Compatible
VWR	UnoCycler	96 x 0.2ml	✓	Compatible
VWR	UnoCycler fast	96 x 0.2ml	✓	Compatible
Whatman Biometra	T1 Thermocycler	96 x 0.2ml	✓	Compatible
Whatman Biometra	T1 Thermocycler Combi	96-77 96 x 0.2ml	✓	Compatible
Whatman Biometra	TGradient	96 x 0.2ml	✓	Compatible
Whatman Biometra	TProfessional	96 x 0.2ml	✓	Compatible
Whatman Biometra	TProfessional Basic	96 x 0.2ml	✓	Compatible
Whatman Biometra	TProfessional Basic Gradient	96 x 0.2ml	✓	Compatible
Whatman Biometra	TProfessional Standard Gradient	96 x 0.2ml	✓	Compatible
Whatman Biometra	TRobot	96 x 0.2ml	✓	Compatible
Whatman Biometra	UNO II	96 x 0.2ml	✓	Compatible
Whatman Biometra	UNO-thermoblock	96 x 0.2ml	✓	Compatible