



Expell PCR Plates Compatibility Chart
Thermocycler Compatibility for 5100315C
ref. PIP_DS_5100315C_22032019

Agilent Technologies	AriaMx Realtime PCR System	96 x 0.2ml		Compatible
Applied Biosystems	7500 Fast Real-Time PCR System	96 x 0.2ml Fast		Compatible
Applied Biosystems	7900HT Fast Real-Time PCR System	96 x 0.2ml Fast		Compatible
Applied Biosystems	9800 Fast Thermal Cyclers	96 x 0.2ml		Compatible
Applied Biosystems	QuantStudio™ 12k Flex	96 x 0.1ml Fast		Compatible
Applied Biosystems	QuantStudio™ 3 Real-Time PCR System	96 x 0.1ml Fast		Compatible
Applied Biosystems	QuantStudio™ 5 Real-Time PCR System	96 x 0.1ml Fast		Compatible
Applied Biosystems	QuantStudio™ 6 Flex	96 x 0.1ml Fast		Compatible
Applied Biosystems	QuantStudio™ 7 Flex	96 x 0.1ml Fast		Compatible
Applied Biosystems	StepOnePlus™ Real-Time PCR System	96 x 0.1ml		Compatible
Applied Biosystems	Veriti™ 96-Well Fast Thermal Cyclers	96 x 0.1ml		Compatible
Applied Biosystems	Veriti™ Dx 96-Well Fast Thermal Cyclers	96 x 0.1ml		Compatible
Applied Biosystems	ViiA™ 7 Real-Time PCR System Fast	96 x 0.1ml		Compatible
Biomérieux	GENE-UP	96 x 0.2ml Fast		Compatible