



Applied Biosystems	7300 Real-Time PCR System	96 x 0.2ml	compatible
<b>Applied Biosystems</b>	<b>7500 Real-Time PCR System</b>	<b>96 x 0.2ml Standard</b>	<b>Single tubes are not recommended</b>
<b>Applied Biosystems</b>	<b>7900HT Real-Time PCR System</b>	<b>96 x 0.2ml Standard</b>	<b>Tubes &amp; stips do not fit. A grid is required</b>
Applied Biosystems	ABI PRISM 7000 SDS	96 x 0.2ml	compatible
Applied Biosystems	ABI PRISM 7700 SDS	96 x 0.2ml	compatible
Applied Biosystems	GeneAmp PCR System 2400	24 x 0.2ml	compatible
Applied Biosystems	GeneAmp PCR System 2700	96 x 0.2ml	compatible
Applied Biosystems	GeneAmp PCR System 2720	96 x 0.2ml	compatible
Applied Biosystems	GeneAmp PCR System 5700	96 x 0.2ml	compatible
Applied Biosystems	GeneAmp PCR System 9600	96 x 0.2ml	compatible
Applied Biosystems	GeneAmp PCR System 9700	96 x 0.2ml	compatible
Applied Biosystems	GeneAmp PCR System 9700	Dual 96-Well 96 x 0.2ml	compatible
Applied Biosystems	GeneAmp PCR System 9700 Aluminum	96 x 0.2ml	compatible
Applied Biosystems	GeneAmp PCR System 9700 Gold	96 x 0.2ml	compatible
Applied Biosystems	ProFlex™ PCR System	32 x 0.2ml	compatible
Applied Biosystems	ProFlex™ PCR System	96 x 0.2ml	compatible
<b>Applied Biosystems</b>	<b>QuantStudio™ 12k Flex</b>	<b>96 x 0.2ml</b>	<b>Single tubes are not recommended</b>
<b>Applied Biosystems</b>	<b>QuantStudio™ 12k Flex</b>	<b>96 x 0.2ml Fast</b>	<b>Does not fit this thermocycler requires low profile products.</b>
<b>Applied Biosystems</b>	<b>QuantStudio™ 6 Flex</b>	<b>96 x 0.2ml</b>	<b>Single tubes are not recommended</b>
<b>Applied Biosystems</b>	<b>QuantStudio™ 6 Flex</b>	<b>96 x 0.2ml Fast</b>	<b>Does not fit this thermocycler requires low profile products</b>
<b>Applied Biosystems</b>	<b>QuantStudio™ 7 Flex</b>	<b>96 x 0.2ml</b>	<b>Single tubes are not recommended</b>
<b>Applied Biosystems</b>	<b>QuantStudio™ 7 Flex</b>	<b>96 x 0.2ml Fast</b>	<b>Does not fit this thermocycler requires low profile products</b>
Applied Biosystems	SimpliAmp™	96 x 0.2ml	compatible
Applied Biosystems	Solid3	96 x 0.2ml	compatible
Applied Biosystems	Veriti™ 96-Well Thermal Cycler	96 x 0.2ml	compatible
Applied Biosystems	Veriti™ Dx 96-Well Thermal Cycler	96 x 0.2ml	compatible
<b>Applied Biosystems</b>	<b>ViiA™ 7 Real-Time PCR System</b>	<b>96 x 0.2ml</b>	<b>Single tubes are not recommended</b>
Bio-Rad Laboratories	C1000	96-Well Fast	compatible
Bio-Rad Laboratories	C1000	Dual 48/48 Fast	compatible
Bio-Rad Laboratories	C1000 Touch	96-Deep Well	compatible
Bio-Rad Laboratories	CFX96 Real-Time PCR Detection System	96 x 0.2ml	compatible
Bio-Rad Laboratories	DNA Engine	48-48 Alpha unit 48 x 0.2ml	compatible
Bio-Rad Laboratories	DNA Engine	96v Alpha unit 96 x 0.2ml	compatible
Bio-Rad Laboratories	DNA Engine Dyad	48-48 Alpha unit 48 x 0.2ml	compatible
Bio-Rad Laboratories	DNA Engine Dyad	96v Alpha unit 96 x 0.2ml	compatible
Bio-Rad Laboratories	DNA Engine Dyad Disciple	48-48 Alpha unit 48 x 0.2ml	compatible
Bio-Rad Laboratories	DNA Engine Dyad Disciple	96v Alpha unit 96 x 0.2ml	compatible
Bio-Rad Laboratories	DNA Engine Tetrad 2	48-48 Alpha unit 48 x 0.2ml	compatible
Bio-Rad Laboratories	DNA Engine Tetrad 2	96v Alpha unit 96 x 0.2ml	compatible
Bio-Rad Laboratories	GeneCycler	24 x 0.2ml	compatible
Bio-Rad Laboratories	iCycler	48 x 0.2ml	compatible
Bio-Rad Laboratories	iCycler	96 x 0.2ml	compatible
Bio-Rad Laboratories	iQ Real-Time PCR	96 x 0.2ml	compatible
Bio-Rad Laboratories	iQ5 Real-Time PCR	96 x 0.2ml	compatible
Bio-Rad Laboratories	MJ Mini	48-12 48 x 0.2ml	compatible
Bio-Rad Laboratories	MyCycler	96 x 0.2ml	compatible
Bio-Rad Laboratories	MyiQ Real-Time PCR	96 x 0.2ml	compatible
Bio-Rad Laboratories	MyiQ2 Real-Time PCR	96 x 0.2ml	compatible
Bio-Rad Laboratories	PTC-100	96 x 0.2ml	compatible
Bio-Rad Laboratories	S1000	96-Deep Well	compatible
Bio-Rad Laboratories	S1000	96-Well Fast	compatible
Bio-Rad Laboratories	S1000	Dual 48/48 Fast	compatible
Bio-Rad Laboratories	T100	96 x 0.2ml	compatible
Eppendorf	Mastercycler	96 x 0.2ml	compatible
Eppendorf	Mastercycler	96-77 96 x 0.2ml	compatible
Eppendorf	Mastercycler ep gradient	96 x 0.2ml	compatible
Eppendorf	Mastercycler ep gradient Silver	96 x 0.2ml	compatible
Eppendorf	Mastercycler ep realplex	96 x 0.2ml	compatible
Eppendorf	Mastercycler ep realplex Silver	96 x 0.2ml	compatible

Eppendorf	Mastercycler gradient	96 x 0.2ml	compatible
Eppendorf	Mastercycler gradient	96-77 96 x 0.2ml	compatible
Eppendorf	Mastercycler nexus	96 x 0.2ml	compatible
Eppendorf	Mastercycler nexus eco	96 x 0.2ml	compatible
Eppendorf	Mastercycler nexus gradient	96 x 0.2ml	compatible
Eppendorf	Mastercycler nexus gradient eco	96 x 0.2ml	compatible
Eppendorf	Mastercycler nexus GSX1	96 x 0.2ml	compatible
Eppendorf	Mastercycler nexus GSX1e	96 x 0.2ml	compatible
Eppendorf	Mastercycler nexus SX1	96 x 0.2ml	compatible
Eppendorf	Mastercycler nexus SX1e	96 x 0.2ml	compatible
Eppendorf	Mastercycler Personal	25-16 25 x 0.2ml	compatible
Eppendorf	Mastercycler pro	96 x 0.2ml	compatible
Eppendorf	Mastercycler pro S	96 x 0.2ml	compatible
Eppendorf	Thermomixer comfort	96 x 0.2ml	compatible
MJ Research Inc.	PTC-100	96 x 0.2ml	compatible
MJ Research Inc.	PTC-100	96 x 0.2ml AgV	compatible
MJ Research Inc.	PTC-150 MiniCycler	25 x 0.2ml	compatible
MJ Research Inc.	PTC-200 DNA Engine	30-48 Alpha unit 48 x 0.2ml	compatible
MJ Research Inc.	PTC-200 DNA Engine	48-48 Alpha unit 48 x 0.2ml	compatible
MJ Research Inc.	PTC-200 DNA Engine	96v Alpha unit 96 x 0.2ml	compatible
MJ Research Inc.	PTC-200 DNA Engine	96v Moto Alpha Unit 96 x 0.2ml	compatible
MJ Research Inc.	PTC-200 DNA Engine	Chromo4 96 x 0.2ml	compatible
MJ Research Inc.	PTC-220 DNA Engine Dyad	30-48 Alpha unit 48 x 0.2ml	compatible
MJ Research Inc.	PTC-220 DNA Engine Dyad	48-48 Alpha unit 48 x 0.2ml	compatible
MJ Research Inc.	PTC-220 DNA Engine Dyad	96v Alpha unit 96 x 0.2ml	compatible
MJ Research Inc.	PTC-220 DNA Engine Dyad	96v Moto Alpha Unit 96 x 0.2ml	compatible
MJ Research Inc.	PTC-220 DNA Engine Dyad	Chromo4 96 x 0.2ml	compatible
MJ Research Inc.	PTC-221 DNA Engine Disciple	30-48 Alpha unit 48 x 0.2ml	compatible
MJ Research Inc.	PTC-221 DNA Engine Disciple	48-48 Alpha unit 48 x 0.2ml	compatible
MJ Research Inc.	PTC-221 DNA Engine Disciple	96v Alpha unit 96 x 0.2ml	compatible
MJ Research Inc.	PTC-221 DNA Engine Disciple	96v Moto Alpha Unit 96 x 0.2ml	compatible
MJ Research Inc.	PTC-221 DNA Engine Disciple	Chromo4 96 x 0.2ml	compatible
MJ Research Inc.	PTC-225 DNA Engine Tetrad	30-48 Alpha unit 48 x 0.2ml	compatible
MJ Research Inc.	PTC-225 DNA Engine Tetrad	48-48 Alpha unit 48 x 0.2ml	compatible
MJ Research Inc.	PTC-225 DNA Engine Tetrad	96v Alpha unit 96 x 0.2ml	compatible
MJ Research Inc.	PTC-225 DNA Engine Tetrad	96v Moto Alpha Unit 96 x 0.2ml	compatible
MJ Research Inc.	PTC-225 DNA Engine Tetrad	Chromo4 96 x 0.2ml	compatible
MJ Research Inc.	PTC-240 DNA Engine Tetrad 2	30-48 Alpha unit 48 x 0.2ml	compatible
MJ Research Inc.	PTC-240 DNA Engine Tetrad 2	48-48 Alpha unit 48 x 0.2ml	compatible
MJ Research Inc.	PTC-240 DNA Engine Tetrad 2	96v Alpha unit 96 x 0.2ml	compatible
MJ Research Inc.	PTC-240 DNA Engine Tetrad 2	96v Moto Alpha Unit 96 x 0.2ml	compatible
MJ Research Inc.	PTC-240 DNA Engine Tetrad 2	Chromo4 96 x 0.2ml	compatible
Qiagen	Rotor-Gene Q®	Rotor 36	compatible
Stratagene	Mx3000P™	96 x 0.2ml	compatible
Stratagene	Mx3005P™	96 x 0.2ml	compatible
Stratagene	Mx4000™	96 x 0.2ml	compatible
Stratagene	Robocycler Gradient 40	96 x 0.2ml	compatible
Whatman Biometra	T1 Thermocycler	96 x 0.2ml	compatible
Whatman Biometra	T1 Thermocycler Combi	96-77 96 x 0.2ml	compatible
Whatman Biometra	T3 Thermocycler	48 x 0.2ml	compatible
Whatman Biometra	T3 Thermocycler Combi	48-18 48 x 0.2ml	compatible
Whatman Biometra	T3000 Thermocycler	48 x 0.2ml	compatible
Whatman Biometra	T3000 Thermocycler Combi	48-18 48 x 0.2ml	compatible
Whatman Biometra	TGradient	96 x 0.2ml	compatible
Whatman Biometra	TOptical	96 x 0.2ml	compatible
Whatman Biometra	TPersonal	48 x 0.2ml	compatible
Whatman Biometra	TPersonal Combi	48-18 48 x 0.2ml	compatible
Whatman Biometra	TProfessional	96 x 0.2ml	compatible
Whatman Biometra	TProfessional Basic	96 x 0.2ml	compatible
Whatman Biometra	TProfessional Basic Gradient	96 x 0.2ml	compatible
Whatman Biometra	TProfessional TRIO	48 x 0.2ml	compatible
<b>Whatman Biometra</b>	<b>TProfessional TRIO combi</b>	<b>48-30 30 x 0.5ml</b>	<b>Does not fit 0.5 ml cyclers</b>
Whatman Biometra	TProfessional TRIO combi	48-30 48 x 0.2ml	compatible
Whatman Biometra	TRobot	96 x 0.2ml	compatible

Whatman Biometra

UNO II

96 x 0.2ml

compatible

Whatman Biometra

UNO-thermoblock

96 x 0.2ml

compatible