



Expell PCR Tubes Compatibility Chart

Thermocycler Compatibility for 5100300C

ref. PIP_DS_5100300C_14012015

Applied Biosystems	7500 Fast Real-Time PCR System	96 x 0.2ml Fast	Single tubes are not recommended
Applied Biosystems	7900HT Fast Real-Time PCR System	96 x 0.2ml Fast	Tubes and strips do not fit
Applied Biosystems	9800 Fast Thermal Cycler	96 x 0.2ml	compatible
Applied Biosystems	GeneAmp PCR System 9600	96 x 0.2ml	compatible
Applied Biosystems	QuantStudio™ 12k Flex	96 x 0.2ml	Does not fit as this cycler has no height adjustable lid
Applied Biosystems	QuantStudio™ 6 Flex	96 x 0.2ml	Does not fit as this cycler has no height adjustable lid
Applied Biosystems	QuantStudio™ 7 Flex	96 x 0.2ml	Does not fit as this cycler has no height adjustable lid
Applied Biosystems	StepOnePlus™ Real-Time PCR System	96 x 0.1ml	compatible
Applied Biosystems	StepOne™ Real-Time PCR System	48 x 0.1ml	compatible
Applied Biosystems	Veriti™ 96-Well Fast Thermal Cycler	96 x 0.1ml	compatible
Applied Biosystems	Veriti™ Dx 96-Well Fast Thermal Cycler	96 x 0.1ml	compatible
Applied Biosystems	ViiA™ 7 Real-Time PCR System	96 x 0.2ml	Does not fit as this cycler has no height adjustable lid
Bio-Rad Laboratories	C1000	96-Well Fast	compatible
Bio-Rad Laboratories	C1000	Dual 48/48 Fast	compatible
Bio-Rad Laboratories	C1000 Touch	96-Deep Well	compatible
Bio-Rad Laboratories	C1000 Touch	96-Well Fast	compatible
Bio-Rad Laboratories	C1000 Touch	Dual 48/48 Fast	compatible
Bio-Rad Laboratories	CFX96 Real-Time PCR Detection System	96 x 0.2ml	compatible
Bio-Rad Laboratories	DNA Engine	48-48 Alpha unit 48 x 0.2ml	compatible
Bio-Rad Laboratories	DNA Engine	96v Alpha unit 96 x 0.2ml	compatible
Bio-Rad Laboratories	DNA Engine	Chromo4 96 x 0.2ml	compatible
Bio-Rad Laboratories	DNA Engine Dyad	48-48 Alpha unit 48 x 0.2ml	compatible
Bio-Rad Laboratories	DNA Engine Dyad	96v Alpha unit 96 x 0.2ml	compatible
Bio-Rad Laboratories	DNA Engine Dyad	Chromo4 96 x 0.2ml	compatible
Bio-Rad Laboratories	DNA Engine Dyad Disciple	48-48 Alpha unit 48 x 0.2ml	compatible
Bio-Rad Laboratories	DNA Engine Dyad Disciple	96v Alpha unit 96 x 0.2ml	compatible
Bio-Rad Laboratories	DNA Engine Dyad Disciple	Chromo4 96 x 0.2ml	compatible
Bio-Rad Laboratories	DNA Engine Tetrad 2	48-48 Alpha unit 48 x 0.2ml	compatible
Bio-Rad Laboratories	DNA Engine Tetrad 2	96v Alpha unit 96 x 0.2ml	compatible
Bio-Rad Laboratories	DNA Engine Tetrad 2	Chromo4 96 x 0.2ml	compatible
Bio-Rad Laboratories	MiniOpticon	48-12 48 x 0.2ml	compatible
Bio-Rad Laboratories	MJ Mini	48-12 12 x 0.5ml	compatible
Bio-Rad Laboratories	MJ Mini	48-12 48 x 0.2ml	compatible
Bio-Rad Laboratories	Opticon 2	96 x 0.2ml	compatible
Bio-Rad Laboratories	PTC-100	96 x 0.2ml	compatible
Bio-Rad Laboratories	S1000	96-Deep Well	compatible
Bio-Rad Laboratories	S1000	96-Well Fast	compatible
Bio-Rad Laboratories	S1000	Dual 48/48 Fast	compatible
Eppendorf	Mastercycler	96 x 0.2ml	compatible
Eppendorf	Mastercycler	96-77 96 x 0.2ml	compatible
Eppendorf	Mastercycler ep gradient	96 x 0.2ml	compatible
Eppendorf	Mastercycler ep gradient Silver	96 x 0.2ml	compatible
Eppendorf	Mastercycler ep realplex	96 x 0.2ml	compatible
Eppendorf	Mastercycler ep realplex Silver	96 x 0.2ml	compatible
Eppendorf	Mastercycler gradient	96 x 0.2ml	compatible
Eppendorf	Mastercycler gradient	96-77 96 x 0.2ml	compatible
Eppendorf	Mastercycler nexus	96 x 0.2ml	compatible

Eppendorf	Mastercycler nexus eco	96 x 0.2ml	compatible
Eppendorf	Mastercycler nexus gradient	96 x 0.2ml	compatible
Eppendorf	Mastercycler nexus gradient eco	96 x 0.2ml	compatible
Eppendorf	Mastercycler Personal	25-16 25 x 0.2ml	compatible
Eppendorf	Mastercycler pro	96 x 0.2ml	compatible
Eppendorf	Mastercycler pro S	96 x 0.2ml	compatible
Eppendorf	Thermomixer comfort	96 x 0.2ml	compatible
MJ Research Inc.	Opticon 2 System	96 x 0.2ml	compatible
MJ Research Inc.	Opticon System	96 x 0.2ml	compatible
MJ Research Inc.	PTC-100	96 x 0.2ml	compatible
MJ Research Inc.	PTC-100	96 x 0.2ml AgV	compatible
MJ Research Inc.	PTC-150 MiniCycler	25 x 0.2ml	compatible
MJ Research Inc.	PTC-200 DNA Engine	30-48 Alpha unit 48 x 0.2ml	compatible
MJ Research Inc.	PTC-200 DNA Engine	48-48 Alpha unit 48 x 0.2ml	compatible
MJ Research Inc.	PTC-200 DNA Engine	96v Alpha unit 96 x 0.2ml	compatible
MJ Research Inc.	PTC-200 DNA Engine	96v Moto Alpha Unit 96 x 0.2ml	compatible
MJ Research Inc.	PTC-200 DNA Engine	Chromo4 96 x 0.2ml	compatible
MJ Research Inc.	PTC-220 DNA Engine Dyad	30-48 Alpha unit 48 x 0.2ml	compatible
MJ Research Inc.	PTC-220 DNA Engine Dyad	48-48 Alpha unit 48 x 0.2ml	compatible
MJ Research Inc.	PTC-220 DNA Engine Dyad	96v Alpha unit 96 x 0.2ml	compatible
MJ Research Inc.	PTC-220 DNA Engine Dyad	96v Moto Alpha Unit 96 x 0.2ml	compatible
MJ Research Inc.	PTC-220 DNA Engine Dyad	Chromo4 96 x 0.2ml	compatible
MJ Research Inc.	PTC-221 DNA Engine Disciple	30-48 Alpha unit 48 x 0.2ml	compatible
MJ Research Inc.	PTC-221 DNA Engine Disciple	48-48 Alpha unit 48 x 0.2ml	compatible
MJ Research Inc.	PTC-221 DNA Engine Disciple	96v Alpha unit 96 x 0.2ml	compatible
MJ Research Inc.	PTC-221 DNA Engine Disciple	96v Moto Alpha Unit 96 x 0.2ml	compatible
MJ Research Inc.	PTC-221 DNA Engine Disciple	Chromo4 96 x 0.2ml	compatible
MJ Research Inc.	PTC-225 DNA Engine Tetrad	30-48 Alpha unit 48 x 0.2ml	compatible
MJ Research Inc.	PTC-225 DNA Engine Tetrad	48-48 Alpha unit 48 x 0.2ml	compatible
MJ Research Inc.	PTC-225 DNA Engine Tetrad	96v Alpha unit 96 x 0.2ml	compatible
MJ Research Inc.	PTC-225 DNA Engine Tetrad	96v Moto Alpha Unit 96 x 0.2ml	compatible
MJ Research Inc.	PTC-225 DNA Engine Tetrad	Chromo4 96 x 0.2ml	compatible
MJ Research Inc.	PTC-240 DNA Engine Tetrad 2	30-48 Alpha unit 48 x 0.2ml	compatible
MJ Research Inc.	PTC-240 DNA Engine Tetrad 2	48-48 Alpha unit 48 x 0.2ml	compatible
MJ Research Inc.	PTC-240 DNA Engine Tetrad 2	96v Alpha unit 96 x 0.2ml	compatible
MJ Research Inc.	PTC-240 DNA Engine Tetrad 2	96v Moto Alpha Unit 96 x 0.2ml	compatible
MJ Research Inc.	PTC-240 DNA Engine Tetrad 2	Chromo4 96 x 0.2ml	compatible
Roche	LightCycler @ 480 Real-Time PCR System	96 x 0.2ml	Do not fit
Roche	LightCycler @ 96 Real-Time PCR System	96 x 0.2ml	compatible
Roche	LightCycler @ Nano	32 x 0.2ml	compatible
Whatman Biometra	T1 Thermocycler	96 x 0.2ml	compatible
Whatman Biometra	T1 Thermocycler Combi	96-77 96 x 0.2ml	compatible
Whatman Biometra	T3 Thermocycler	48 x 0.2ml	compatible
Whatman Biometra	T3 Thermocycler Combi	48-18 48 x 0.2ml	compatible
Whatman Biometra	T3000 Thermocycler	48 x 0.2ml	compatible
Whatman Biometra	T3000 Thermocycler Combi	48-18 48 x 0.2ml	compatible
Whatman Biometra	TGradient	96 x 0.2ml	compatible
Whatman Biometra	TPersonal	48 x 0.2ml	compatible
Whatman Biometra	TPersonal Combi	48-18 48 x 0.2ml	compatible
Whatman Biometra	TProfessional	96 x 0.2ml	compatible
Whatman Biometra	TProfessional Basic	96 x 0.2ml	compatible

Whatman Biometra	TProfessional Basic Gradient	96 x 0.2ml	compatible
Whatman Biometra	TRobot	96 x 0.2ml	compatible
Whatman Biometra	UNO II	96 x 0.2ml	compatible
Whatman Biometra	UNO-thermoblock	96 x 0.2ml	compatible