



Expell PCR Tubes Compatibility Chart

Thermocycler Compatibility for 5100311C

ref. PIP_DS_5100311C_14012015

Applied Biosystems	7500 Fast Real-Time PCR System	96 x 0.2ml Fast	compatible
Applied Biosystems	7900HT Fast Real-Time PCR System	96 x 0.2ml Fast	Do not fit
Applied Biosystems	7900HT Real-Time PCR System	384-well	Does not fit as this cyclcr is a 384 well cyclcr
Applied Biosystems	9800 Fast Thermal Cyclcr	96 x 0.2ml	compatible
Applied Biosystems	GeneAmp PCR System 9600	96 x 0.2ml	compatible
Applied Biosystems	QuantStudio™ 12k Flex	384-well	Does not fit as this cyclcr is a 384 well cyclcr
Applied Biosystems	QuantStudio™ 12k Flex	96 x 0.2ml Fast	Single tubes and strips are not recommended
Applied Biosystems	QuantStudio™ 6 Flex	384-well	Does not fit as this cyclcr is a 384 well cyclcr
Applied Biosystems	QuantStudio™ 6 Flex	96 x 0.2ml Fast	Single tubes and strips are not recommended
Applied Biosystems	QuantStudio™ 7 Flex	384-well	Does not fit as this cyclcr is a 384 well cyclcr
Applied Biosystems	QuantStudio™ 7 Flex	96 x 0.2ml Fast	Single tubes and strips are not recommended
Applied Biosystems	StepOnePlus™ Real-Time PCR System	96 x 0.1ml	compatible
Applied Biosystems	StepOne™ Real-Time PCR System	48 x 0.1ml	compatible
Applied Biosystems	Veriti™ 96-Well Fast Thermal Cyclcr	96 x 0.1ml	compatible
Applied Biosystems	Veriti™ Dx 96-Well Fast Thermal Cyclcr	96 x 0.1ml	compatible
Applied Biosystems	ViiA™ 7 Real-Time PCR System Fast	96 x 0.1ml	Single tubes and strips are not recommended
Bio-Rad Laboratories	C1000	96-Well Fast	compatible
Bio-Rad Laboratories	C1000	Dual 48/48 Fast	compatible
Bio-Rad Laboratories	C1000 Touch	96-Deep Well	compatible
Bio-Rad Laboratories	C1000 Touch	96-Well Fast	compatible
Bio-Rad Laboratories	C1000 Touch	Dual 48/48 Fast	compatible
Bio-Rad Laboratories	DNA Engine	48-48 Alpha unit 48 x 0.2ml	compatible
Bio-Rad Laboratories	DNA Engine	96v Alpha unit 96 x 0.2ml	compatible
Bio-Rad Laboratories	DNA Engine Dyad	48-48 Alpha unit 48 x 0.2ml	compatible
Bio-Rad Laboratories	DNA Engine Dyad	96v Alpha unit 96 x 0.2ml	compatible
Bio-Rad Laboratories	DNA Engine Dyad	Chromo4 96 x 0.2ml	compatible
Bio-Rad Laboratories	DNA Engine Dyad Disciple	48-48 Alpha unit 48 x 0.2ml	compatible
Bio-Rad Laboratories	DNA Engine Dyad Disciple	96v Alpha unit 96 x 0.2ml	compatible
Bio-Rad Laboratories	DNA Engine Dyad Disciple	Chromo4 96 x 0.2ml	compatible
Bio-Rad Laboratories	DNA Engine Tetrad 2	48-48 Alpha unit 48 x 0.2ml	compatible
Bio-Rad Laboratories	DNA Engine Tetrad 2	96v Alpha unit 96 x 0.2ml	compatible
Bio-Rad Laboratories	DNA Engine Tetrad 2	Chromo4 96 x 0.2ml	compatible
Bio-Rad Laboratories	MJ Mini	48-12 48 x 0.2ml	compatible
Bio-Rad Laboratories	PTC-100	96 x 0.2ml	compatible
Bio-Rad Laboratories	S1000	96-Deep Well	compatible
Bio-Rad Laboratories	S1000	96-Well Fast	compatible
Bio-Rad Laboratories	S1000	Dual 48/48 Fast	compatible
Eppendorf	Mastercycler	96 x 0.2ml	compatible
Eppendorf	Mastercycler	96-77 96 x 0.2ml	compatible
Eppendorf	Mastercycler ep gradient	96 x 0.2ml	compatible
Eppendorf	Mastercycler ep gradient Silver	96 x 0.2ml	compatible
Eppendorf	Mastercycler gradient	96 x 0.2ml	compatible
Eppendorf	Mastercycler gradient	96-77 96 x 0.2ml	compatible
Eppendorf	Mastercycler nexus	96 x 0.2ml	compatible

Eppendorf	Mastercycler nexus eco	96 x 0.2ml	compatible
Eppendorf	Mastercycler nexus gradient	96 x 0.2ml	compatible
Eppendorf	Mastercycler nexus gradient eco	96 x 0.2ml	compatible
Eppendorf	Mastercycler pro	96 x 0.2ml	compatible
Eppendorf	Mastercycler pro S	96 x 0.2ml	compatible
MJ Research Inc.	PTC-100	96 x 0.2ml	compatible
MJ Research Inc.	PTC-100	96 x 0.2ml AgV	compatible
MJ Research Inc.	PTC-200 DNA Engine	30-48 Alpha unit 30 x 0.5ml	Does not fit in a 0.5 ml block
MJ Research Inc.	PTC-200 DNA Engine	48-48 Alpha unit 48 x 0.2ml	compatible
MJ Research Inc.	PTC-200 DNA Engine	96v Alpha unit 96 x 0.2ml	compatible
MJ Research Inc.	PTC-200 DNA Engine	96v Moto Alpha Unit 96 x 0.2ml	compatible
MJ Research Inc.	PTC-220 DNA Engine Dyad	30-48 Alpha unit 48 x 0.2ml	compatible
MJ Research Inc.	PTC-220 DNA Engine Dyad	48-48 Alpha unit 48 x 0.2ml	compatible
MJ Research Inc.	PTC-220 DNA Engine Dyad	96v Alpha unit 96 x 0.2ml	compatible
MJ Research Inc.	PTC-220 DNA Engine Dyad	96v Moto Alpha Unit 96 x 0.2ml	compatible
MJ Research Inc.	PTC-221 DNA Engine Disciple	30-48 Alpha unit 48 x 0.2ml	compatible
MJ Research Inc.	PTC-221 DNA Engine Disciple	48-48 Alpha unit 48 x 0.2ml	compatible
MJ Research Inc.	PTC-221 DNA Engine Disciple	96v Alpha unit 96 x 0.2ml	compatible
MJ Research Inc.	PTC-221 DNA Engine Disciple	96v Moto Alpha Unit 96 x 0.2ml	compatible
MJ Research Inc.	PTC-225 DNA Engine Tetrad	30-48 Alpha unit 48 x 0.2ml	compatible
MJ Research Inc.	PTC-225 DNA Engine Tetrad	48-48 Alpha unit 48 x 0.2ml	compatible
MJ Research Inc.	PTC-225 DNA Engine Tetrad	96v Alpha unit 96 x 0.2ml	compatible
MJ Research Inc.	PTC-225 DNA Engine Tetrad	96v Moto Alpha Unit 96 x 0.2ml	compatible
MJ Research Inc.	PTC-240 DNA Engine Tetrad 2	30-48 Alpha unit 48 x 0.2ml	compatible
MJ Research Inc.	PTC-240 DNA Engine Tetrad 2	48-48 Alpha unit 48 x 0.2ml	compatible
MJ Research Inc.	PTC-240 DNA Engine Tetrad 2	96v Alpha unit 96 x 0.2ml	compatible
MJ Research Inc.	PTC-240 DNA Engine Tetrad 2	96v Moto Alpha Unit 96 x 0.2ml	compatible
Roche	LightCycler @ Nano	32 x 0.2ml	compatible
Whatman Biometra	T1 Thermocycler	96 x 0.2ml	compatible
Whatman Biometra	T1 Thermocycler Combi	96-77 96 x 0.2ml	compatible
Whatman Biometra	T3 Thermocycler	48 x 0.2ml	compatible
Whatman Biometra	T3 Thermocycler Combi	48-18 48 x 0.2ml	compatible
Whatman Biometra	T3000 Thermocycler	48 x 0.2ml	compatible
Whatman Biometra	T3000 Thermocycler Combi	48-18 48 x 0.2ml	compatible
Whatman Biometra	TGradient	96 x 0.2ml	compatible
Whatman Biometra	TOptical	96 x 0.2ml	compatible
Whatman Biometra	TPersonal	48 x 0.2ml	compatible
Whatman Biometra	TPersonal Combi	48-18 48 x 0.2ml	compatible
Whatman Biometra	TProfessional	96 x 0.2ml	compatible
Whatman Biometra	TProfessional Basic	96 x 0.2ml	compatible
Whatman Biometra	TProfessional Basic Gradient	96 x 0.2ml	compatible
Whatman Biometra	TProfessional TRIO	48 x 0.2ml	compatible
Whatman Biometra	TProfessional TRIO combi	48-30 48 x 0.2ml	compatible
Whatman Biometra	TRobot	96 x 0.2ml	compatible
Whatman Biometra	UNO II	96 x 0.2ml	compatible
Whatman Biometra	UNO-thermoblock	96 x 0.2ml	compatible