

COMPATIBILITY CHART

| | | CAPPSolo | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
|---|----------|--|--------------|--------------|--------------|--------------|----------------|----------------|----------------|----------------|------------------|---|-----------------|---------------|-----------------|-----------------|-----------------|---------------|---------------|---|------------------|---------------------|--|------------------|-----------------|------------------|------------------|------------------|-------------------|-------------------|---|
| | | CAPPSolo single channel pipette fixed volume | | | | | | | | | | CAPPSolo single channel pipette variable volume | | | | | | | | CAPPSolo multichannel pipette 8-channel | | | CAPPSolo multichannel pipette 12-channel | | | | | | | | |
| | | 5µl/\$5-1F | 10µl/\$10-1F | 20µl/\$20-1F | 25µl/\$25-1F | 50µl/\$50-1F | 100µl/\$100-1F | 200µl/\$200-1F | 250µl/\$250-1F | 500µl/\$500-1F | 1000µl/\$1000-1F | 2ml/\$2000-1F | 2.5ml/\$2500-1F | 5ml/\$5000-1F | 10ml/\$10000-1F | 0.1-2.5µl/\$2-1 | 0.5-10µl/\$10-1 | 2-20µl/\$20-1 | 5-50µl/\$50-1 | 10-100µl/\$100-1 | 20-200µl/\$200-1 | 100-1000µl/\$1000-1 | 0.5-5ml/\$5000-1 | 2-10ml/\$10000-1 | 0.5-10µl/\$10-8 | 10-100µl/\$100-8 | 30-300µl/\$300-8 | 0.5-10µl/\$10-12 | 10-100µl/\$100-12 | 30-300µl/\$300-12 | |
| Description | Cat. No. | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| CAPP Expell Filter Tips | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Expell 10 µL, w/ filter, clear, sterile | 5130030C | | | | | | | | | | | | | | ■ | ■ | ■ | | | | | | | | | ■ | | | ■ | | |
| Expell 10 µL, w/ filter, clear | 5131030C | | | | | | | | | | | | | | ■ | ■ | ■ | | | | | | | | | ■ | | | ■ | | |
| Expell 10 µL XL, w/ filter, clear, sterile | 5130060C | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Expell 10 µL XL, w/ filter, clear | 5131060C | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Expell 10/20 µL XL, extra narrow, w/ filter, clear, sterile | 5130061C | | | | | | | | | | | | | | ■ | ■ | ■ | | | | | | | | | ■ | | | ■ | | |
| Expell 10/20 µL XL, extra narrow, w/ filter, clear | 5131061C | | | | | | | | | | | | | | ■ | ■ | ■ | | | | | | | | ■ | | | ■ | | | |
| Expell 20 µL, w/ filter, clear, sterile | 5130062C | ■ | ■ | ■ | ■ | ■ | ■ | ■ | | | | | | | | | | ■ | ■ | ■ | | | | | | | ■ | ■ | | ■ | ■ |
| Expell 20 µL, w/ filter, clear | 5131062C | ■ | ■ | ■ | ■ | ■ | ■ | ■ | | | | | | | | | | ■ | ■ | ■ | | | | | | | ■ | ■ | | ■ | ■ |
| Expell 100 µL, w/ filter, clear, sterile | 5130066C | ■ | ■ | ■ | ■ | ■ | ■ | ■ | | | | | | | | | | ■ | ■ | ■ | | | | | | | ■ | ■ | | ■ | ■ |
| Expell 100 µL, w/ filter, clear | 5131066C | ■ | ■ | ■ | ■ | ■ | ■ | ■ | | | | | | | | | | ■ | ■ | ■ | | | | | | | ■ | ■ | | ■ | ■ |
| Expell 200 µL, w/ filter, clear, sterile | 5130090C | | | | | | ■ | ■ | | | | | | | | | | | | ■ | ■ | | | | | | ■ | ■ | | ■ | ■ |
| Expell 200 µL, w/ filter, clear | 5131090C | | | | | | ■ | ■ | | | | | | | | | | | | ■ | ■ | | | | | | ■ | ■ | | ■ | ■ |
| Expell 1000 µL, w/ filter, clear, sterile | 5130150C | | | | | | | | ■ | ■ | ■ | | | | | | | | | | | | ■ | | | | | | | | |
| Expell 1000 µL, w/ filter, clear | 5131150C | | | | | | | | ■ | ■ | ■ | | | | | | | | | | | | ■ | | | | | | | | |
| Expell 1250 µL XL, w/ filter, clear, sterile | 5130123C | | | | | | | | ■ | ■ | ■ | | | | | | | | | | | | ■ | | | | | | | | |
| CAPP Expell and ExpellPlus individually wrapped Pipette Tips | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| ExpellPlus 10 µL, clear, sterile, individually wrapped | 5030012C | | | | | | | | | | | | | | ■ | ■ | ■ | | | | | | | | | ■ | | | ■ | | |
| ExpellPlus 10 µL XL, clear, sterile, individually wrapped | 5030042C | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| ExpellPlus 200 µL, clear, sterile, individually wrapped | 5030072C | ■ | ■ | ■ | ■ | ■ | ■ | ■ | | | | | | | | | | ■ | ■ | ■ | | | | | | | ■ | ■ | | ■ | ■ |
| Expell 1000 µL (1250 µL), clear, sterile, individually wrapped | 5130132C | | | | | | | | ■ | ■ | ■ | | | | | | | | | | | | | | | | | | | | |
| CAPP Expell and ExpellPlus Filter Tips in PaperBox Racks | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| ExpellPlus 10 µL, w/ filter, clear, sterile | 5030029C | | | | | | | | | | | | | | ■ | ■ | ■ | | | | | | | | | ■ | | | ■ | | |
| ExpellPlus 10/20 µL XL, extra narrow, w/ filter, clear, sterile | 5030059C | | | | | | | | | | | | | | ■ | ■ | ■ | | | | | | | | | ■ | | | ■ | | |
| ExpellPlus 20 µL, w/ filter, clear, sterile | 5030063C | ■ | ■ | ■ | ■ | ■ | ■ | ■ | | | | | | | | | | ■ | ■ | ■ | | | | | | | ■ | ■ | | ■ | ■ |
| ExpellPlus 100 µL, w/ filter, clear, sterile | 5030064C | ■ | ■ | ■ | ■ | ■ | ■ | ■ | | | | | | | | | | ■ | ■ | ■ | | | | | | | ■ | ■ | | ■ | ■ |
| ExpellPlus 200 µL, w/ filter, clear, sterile | 5030089C | | | | | | ■ | ■ | | | | | | | | | | | | ■ | ■ | | | | | | ■ | ■ | | ■ | ■ |
| ExpellPlus 300 µL, w/ filter, clear, sterile | 5030114C | | | | | | ■ | ■ | | | | | | | | | | | | ■ | ■ | | | | | | ■ | ■ | | ■ | ■ |
| Expell 1000 µL, w/ filter, clear, sterile | 5130149C | | | | | | | | ■ | ■ | ■ | | | | | | | | | | | | ■ | | | | | | | | |
| ExpellPlus 1250 µL XL, w/ filter, clear, sterile | 5130119C | | | | | | | | ■ | ■ | ■ | | | | | | | | | | | | ■ | | | | | | | | |
| CAPP Expell Pipette Tips Reload System | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Expell 10 µL, clear, reload | 5030035C | | | | | | | | | | | | | | ■ | ■ | ■ | | | | | | | | | ■ | | | ■ | | |
| Expell 10 µL XL, extra narrow, clear, reload | 5030055C | | | | | | | | | | | | | | ■ | ■ | ■ | | | | | | | | | ■ | | | ■ | | |
| Expell 200 µL, clear, reload | 5030095C | ■ | ■ | ■ | ■ | ■ | ■ | ■ | | | | | | | | | | ■ | ■ | ■ | | | | | | | ■ | ■ | | ■ | ■ |
| Expell 300 µL, clear, reload | 5130111C | ■ | ■ | ■ | ■ | ■ | ■ | ■ | | | | | | | | | | ■ | ■ | ■ | | | | | | | ■ | ■ | | ■ | ■ |
| Expell 1000 µL (1250 µL), clear, reload | 5130155C | | | | | | | | ■ | ■ | ■ | | | | | | | | | | | | ■ | | | | | | | | |

■ compatible x not compatible





COMPATIBILITY CHART

| Description | Cat. No. | CAPPSolo | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
|---|----------|--|-------------|-------------|-------------|-------------|---------------|---------------|---------------|---------------|-----------------|--------------|----------------|---|----------------|----------------|----------------|--------------|--------------|---|-----------------|--------------------|--|-----------------|----------------|-----------------|-----------------|-----------------|------------------|------------------|---|
| | | CAPPSolo single channel pipette fixed volume | | | | | | | | | | | | CAPPSolo single channel pipette variable volume | | | | | | CAPPSolo multichannel pipette 8-channel | | | CAPPSolo multichannel pipette 12-channel | | | | | | | | |
| | | 5µL/S5-1F | 10µL/S10-1F | 20µL/S20-1F | 25µL/S25-1F | 50µL/S50-1F | 100µL/S100-1F | 200µL/S200-1F | 250µL/S250-1F | 500µL/S500-1F | 1000µL/S1000-1F | 2ml/S2000-1F | 2.5ml/S2500-1F | 5ml/S5000-1F | 10ml/S10000-1F | 0.1-2.5µL/S2-1 | 0.5-10µL/S10-1 | 2-20µL/S20-1 | 5-50µL/S50-1 | 10-100µL/S100-1 | 20-200µL/S200-1 | 100-1000µL/S1000-1 | 0.5-5ml/S5000-1 | 2-10ml/S10000-1 | 0.5-10µL/S10-8 | 10-100µL/S100-8 | 30-300µL/S300-8 | 0.5-10µL/S10-12 | 10-100µL/S100-12 | 30-300µL/S300-12 | |
| CAPP ExpellPlus Low Retention Pipette Tips | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| ExpellPlus 10 µL, clear, sterile | 5030010C | | | | | | | | | | | | | | | ■ | ■ | ■ | | | | | | | | ■ | | | ■ | | |
| ExpellPlus 10 µL, clear, sterile | 5030020C | | | | | | | | | | | | | | | ■ | ■ | ■ | | | | | | | | ■ | | | ■ | | |
| ExpellPlus 10 µL XL, clear, sterile | 5030040C | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| ExpellPlus 10 µL XL, clear, sterile | 5030050C | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| ExpellPlus 10/20 µL XL, extra narrow, clear, sterile | 5030051C | | | | | | | | | | | | | | | ■ | ■ | ■ | | | | | | | | ■ | | | ■ | | |
| ExpellPlus 200 µL, clear, sterile | 5030070C | ■ | ■ | ■ | ■ | ■ | ■ | ■ | | | | | | | | | | | ■ | ■ | ■ | | | | | | ■ | ■ | | ■ | ■ |
| ExpellPlus 200 µL, clear, sterile | 5030080C | ■ | ■ | ■ | ■ | ■ | ■ | ■ | | | | | | | | | | | ■ | ■ | ■ | | | | | | ■ | ■ | | ■ | ■ |
| ExpellPlus 300 µL, clear, sterile | 5130110C | ■ | ■ | ■ | ■ | ■ | ■ | ■ | | | | | | | | | | | ■ | ■ | ■ | | | | | | ■ | ■ | | ■ | ■ |
| CAPP ExpellPlus Low Retention Filter Tips | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| ExpellPlus 10 µL, w/ filter, clear, sterile | 5030030C | | | | | | | | | | | | | | | ■ | ■ | ■ | | | | | | | | ■ | | | ■ | | |
| ExpellPlus 10 µL XL, w/ filter, clear, sterile | 5030060C | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| ExpellPlus 10/20 µL XL, extra narrow, w/ filter, clear, sterile | 5030061C | | | | | | | | | | | | | | | ■ | ■ | ■ | | | | | | | | ■ | | | ■ | | |
| ExpellPlus 20 µL, w/ filter, clear, sterile | 5030062C | ■ | ■ | ■ | ■ | ■ | ■ | ■ | | | | | | | | | | | ■ | ■ | ■ | | | | | | ■ | ■ | | ■ | ■ |
| ExpellPlus 100 µL, w/ filter, clear, sterile | 5030066C | ■ | ■ | ■ | ■ | ■ | ■ | ■ | | | | | | | | | | | ■ | ■ | ■ | | | | | | ■ | ■ | | ■ | ■ |
| ExpellPlus 200 µL, w/ filter, clear, sterile | 5030090C | | | | | | | | | | | | | | | | | | ■ | ■ | ■ | | | | | | ■ | ■ | | ■ | ■ |
| ExpellPlus 300 µL, w/ filter, clear, sterile | 5030115C | | | | | | | | | | | | | | | | | | ■ | ■ | ■ | | | | | | ■ | ■ | | ■ | ■ |
| ExpellPlus 1000 µL, w/ filter, clear, sterile | 5030150C | | | | | | | | | | ■ | ■ | ■ | | | | | | | | | ■ | | | | | | | | | |
| CAPP Expell Pipette Tips | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Expell 10 µL, clear | 5030015C | | | | | | | | | | | | | | | ■ | ■ | ■ | | | | | | | | ■ | | | ■ | | |
| Expell 10 µL, clear | 5130010C | | | | | | | | | | | | | | | ■ | ■ | ■ | | | | | | | | ■ | | | ■ | | |
| Expell 10 µL XL, clear | 5030045C | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Expell 10 µL XL, clear | 5130040C | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Expell 10/20 µL XL, extra narrow, clear | 5030046C | | | | | | | | | | | | | | | ■ | ■ | ■ | | | | | | | | ■ | | | ■ | | |
| Expell 10/20 µL XL, extra narrow, clear | 5130041C | | | | | | | | | | | | | | | ■ | ■ | ■ | | | | | | | ■ | | | ■ | | | |
| Expell 200 µL, clear | 5030075C | ■ | ■ | ■ | ■ | ■ | ■ | ■ | | | | | | | | | | | ■ | ■ | ■ | | | | | | ■ | ■ | | ■ | ■ |
| Expell 200 µL, clear | 5130070C | ■ | ■ | ■ | ■ | ■ | ■ | ■ | | | | | | | | | | | ■ | ■ | ■ | | | | | | ■ | ■ | | ■ | ■ |
| Expell 200 µL, yellow | 4130075C | | | | | | | | | | | | | | | | | | ■ | ■ | ■ | | | | | | ■ | ■ | | ■ | ■ |
| Expell 300 µL, clear | 5130105C | ■ | ■ | ■ | ■ | ■ | ■ | ■ | | | | | | | | | | | ■ | ■ | ■ | | | | | | ■ | ■ | | ■ | ■ |
| Expell 300 µL, clear | 5130100C | ■ | ■ | ■ | ■ | ■ | ■ | ■ | | | | | | | | | | | ■ | ■ | ■ | | | | | | ■ | ■ | | ■ | ■ |
| Expell 1000 µL (1250 µL), clear | 5130135C | | | | | | | | | | ■ | ■ | ■ | | | | | | | | | ■ | | | | | | | | | |
| Expell 1000 µL (1250 µL), clear | 5130130C | | | | | | | | | | ■ | ■ | ■ | | | | | | | | | ■ | | | | | | | | | |
| Expell 1000 µL (1250 µL), blue | 4130135C | | | | | | | | | | ■ | ■ | ■ | | | | | | | | | ■ | | | | | | | | | |
| Expell 1000 µL (1250 µL), clear, sterile | 5130140C | | | | | | | | | | ■ | ■ | ■ | | | | | | | | | ■ | | | | | | | | | |
| Expell 1250 µL XL, clear | 5130121C | | | | | | | | | | ■ | ■ | ■ | | | | | | | | | ■ | | | | | | | | | |
| Expell 1250 µL XL, clear | 5130120C | | | | | | | | | | ■ | ■ | ■ | | | | | | | | | ■ | | | | | | | | | |
| Expell 1250 µL XL, clear, sterile | 5130122C | | | | | | | | | | ■ | ■ | ■ | | | | | | | | | ■ | | | | | | | | | |
| Expell 5000 µL, clear, sterile | 5130165C | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Expell 5000 µL, clear | 5130160C | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Expell 5000 µL, clear | 5130161C | | | | | | | | | | | ■ | ■ | ■ | | | | | | | | | | | ■ | | | | | | |
| Expell 5500 µL, clear | 5130192C | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Expell 10000 µL, clear, sterile | 5130195C | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Expell 10000 µL, clear | 5130190C | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Expell 10000 µL, clear | 5130191C | | | | | | | | | | | | | | | ■ | | | | | | | | | | ■ | | | | | |

■ compatible x not compatible